


SPECIALTY INTUBATION

SHERIDAN® LTS™ MICROLARYNGEAL ETT

Description:

The Sheridan® LTS™ Microlaryngeal Endotracheal Tube (ETT) features:

- Small diameter (4, 5 or 6 mm) provides better visualization and access to surgical site
- Cuff diameter and tube length is that of standard 8 mm tube
- Designed for microlaryngeal or tracheal surgery
- May also be useful for patients whose airways have been narrowed, preventing the use of a larger tube
- Single-use, sterile



United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699

Lat
Am

Catalog #	Size	Qty/Box
5-11108	Sheridan® LTS™ Microlaryngeal ETT, Size 4.0	10
5-11110	Sheridan® LTS™ Microlaryngeal ETT, Size 5.0	10
5-11112	Sheridan® LTS™ Microlaryngeal ETT, Size 6.0	10

© 2026 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.