


ENDOBRONCHIAL TUBES

SHERIDAN® SHER-I-BRONCH® LEFT-SIDED ENDOBRONCHIAL TUBE

Description:

The Sheridan® Sher-I-Bronch® Left-Sided Endobronchial Tube features:

- Double-lumen endobronchial tubes designed for selective one-lung ventilation, allowing for quiet surgical field on opposite side
- Right-sided tube's double cuff facilitates positioning of tube in right main stem bronchus
- Double-thick tracheal cuff helps to reduce the risk of inadvertent tearing upon insertion
- Single-use, sterile



United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699

Lat

Am

© 2026 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.