


ENDOBRONCHIAL TUBES

SHERIDAN® SHER-I-BRONCH® RIGHT-SIDED ENDOBRONCHIAL TUBE

Description:

The Sheridan® Sher-I-Bronch® Right-Sided Endobronchial Tube features:

- Double-lumen endobronchial tubes designed for selective one-lung ventilation, allowing for quiet surgical field on opposite side
- Right-sided tube’s double cuff facilitates positioning of tube in right main stem bronchus
- Double-thick tracheal cuff helps to reduce the risk of inadvertent tearing upon insertion
- Single-use, sterile



United States

To order this product, please call:
866-246-6990

Canada

To order this product, please call:
800-387-9699

Lat

Am