


ENDOTRACHEAL TUBE ACCESSORIES

SHERIDAN® SHER-I-SWIV/FO™ DOUBLE SWIVEL CONNECTOR

Description:

The Sheridan® Sher-I-SWIV/FO™ Double Swivel Connector for use with fiber optics features:

- Special orifice allows easy passage of fiber optic scope
- Maintains seal around scope during bronchoscopy procedure
- Swivel feature reduces breathing circuit torque on the tracheal tube
- Designed to facilitate introduction of fiber optic bronchoscope during mechanical ventilation



United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699

Lat
Am

Catalog #	Description	Qty/Box
5-15401	Sheridan® Sher-I-SWIV/FO™ Double Swivel Connector	25

© 2026 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.