


SPECIALTY INTUBATION

SHERIDAN® SPIRAL-FLEX® CUFFED ETT

Description:

The Sheridan® Spiral-Flex® Cuffed Endotracheal Tube (ETT) features:

- Spiral steel reinforcing wire in tube wall reduces risk of kinking or compression
- Available with a high-volume, tapered, low pressure, blue cuff
- Cuffed oral tube includes preloaded stylet
- Single-use, sterile



United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699

Lat

Am

Catalog #	Size	Qty/Box
5-12610	Sheridan® Spiral-Flex® Cuffed Oral/Nasal ETT, Size 5.0	5
5-12611	Sheridan® Spiral-Flex® Cuffed Oral/Nasal ETT, Size 5.5	5
5-12612	Sheridan® Spiral-Flex® Cuffed Oral/Nasal ETT, Size 6.0	5
5-12613	Sheridan® Spiral-Flex® Cuffed Oral/Nasal ETT, Size 6.5	5
5-12614	Sheridan® Spiral-Flex® Cuffed Oral/Nasal ETT, Size 7.0	5
5-12615	Sheridan® Spiral-Flex® Cuffed Oral/Nasal ETT, Size 7.5	5

Catalog #	Size	Qty/Box
5-12616	Sheridan® Spiral-Flex® Cuffed Oral/Nasal ETT, Size 8.0	5
5-12617	Sheridan® Spiral-Flex® Cuffed Oral/Nasal ETT, Size 8.5	5
5-12618	Sheridan® Spiral-Flex® Cuffed Oral/Nasal ETT, Size 9.0	5

© 2026 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.