


**SPECIALTY INTUBATION**

**SHERIDAN® SPIRAL-FLEX® MURPHY UNCUFFED ETT**

**Description:**

The Sheridan® Spiral-Flex® Murphy Uncuffed Endotracheal Tube (ETT) features:

- Spiral steel reinforcing wire in tube wall reduces risk of kinking or compression
- Available with a high-volume, tapered, low pressure, blue cuff
- Cuffed oral tube includes preloaded stylet
- Single-use, sterile



United States

To order this product, please call:

**866-246-6990**

Canada

To order this product, please call:

**800-387-9699**

Lat

Am

Catalog #	Size	Qty/Box
5-12906	Sheridan® Spiral-Flex® Murphy Uncuffed ETT, Size 3.0	5
5-12907	Sheridan® Spiral-Flex® Murphy Uncuffed ETT, Size 3.5	5
5-12908	Sheridan® Spiral-Flex® Murphy Uncuffed ETT, Size 4.0	5
5-12909	Sheridan® Spiral-Flex® Murphy Uncuffed ETT, Size 4.5	5
5-12910	Sheridan® Spiral-Flex® Murphy Uncuffed ETT, Size 5.0	5