


SPECIALTY INTUBATION

SHERIDAN® SPIRAL-FLEX® ORAL ETT

Description:

The Sheridan® Spiral-Flex® Oral Endotracheal Tube (ETT) features:

- Spiral steel reinforcing wire in tube wall reduces risk of kinking or compression
- Available with a high-volume, tapered, low pressure, blue cuff
- Cuffed oral tube includes preloaded stylet
- Single-use, sterile



United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699

Lat

Am

Catalog #	Size	Qty/Box
5-12510	Sheridan® Spiral-Flex® Oral ETT, Size 5.0	5
5-12511	Sheridan® Spiral-Flex® Oral ETT, Size 5.5	5
5-12512	Sheridan® Spiral-Flex® Oral ETT, Size 6.0	5
5-12513	Sheridan® Spiral-Flex® Oral ETT, Size 6.5	5
5-12514	Sheridan® Spiral-Flex® Oral ETT, Size 7.0	5
5-12515	Sheridan® Spiral-Flex® Oral ETT, Size 7.5	5

Catalog #	Size	Qty/Box
5-12516	Sheridan® Spiral-Flex® Oral ETT, Size 8.0	5
5-12517	Sheridan® Spiral-Flex® Oral ETT, Size 8.5	5
5-12518	Sheridan® Spiral-Flex® Oral ETT, Size 9.0	5

© 2026 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.