


STANDARD INTUBATION

SHERIDAN UNCUFFED™ MURPHY EYE ETT

Description:

The Sheridan Uncuffed™ Murphy Eye Endotracheal Tube (ETT) is a single-use ETT which softens at body temperature to conform to airway anatomy, and features:

- Smooth, molded, gently cupped tip eases insertion during intubation
- Special depth lines at distal tip reference depth of intubation
- For oral or nasal intubation



United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699

Lat

Am

Catalog #	Size	Qty/Box
5-10404	Sheridan Uncuffed™ Murphy Eye ETT, Size 2.0	10
5-10405	Sheridan Uncuffed™ Murphy Eye ETT, Size 2.5	10
5-10406	Sheridan Uncuffed™ Murphy Eye ETT, Size 3.0	10
5-10407	Sheridan Uncuffed™ Murphy Eye ETT, Size 3.5	10
5-10408	Sheridan Uncuffed™ Murphy Eye ETT, Size 4.0	10
5-10409	Sheridan Uncuffed™ Murphy Eye ETT, Size 4.5	10

Catalog #	Size	Qty/Box
5-10410	Sheridan Uncuffed™ Murphy Eye ETT, Size 5.0	10
5-10411	Sheridan Uncuffed™ Murphy Eye ETT, Size 5.5	10
5-10412	Sheridan Uncuffed™ Murphy Eye ETT, Size 6.0	10
5-10413	Sheridan Uncuffed™ Murphy Eye ETT, Size 6.5	10
5-10414	Sheridan Uncuffed™ Murphy Eye ETT, Size 7.0	10

© 2026 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.