


ETCO₂

SOFTECH® BI-FLO ETCO₂ SAMPLING CANNULA

Description:

The Softech® Bi-Flo ETCO₂ Sampling Cannula allows delivery of supplemental oxygen while simultaneously sampling expired CO₂, featuring a design that allows delivery of oxygen to both nasal prongs and sampling of expired gases from both nares.



United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699

Lat

Am

Catalog #	Size	Tubing Length	Luer Connector	Qty/Case
RHET43U	Softech® Bi-Flo ETCO ₂ Sampling Cannula, Adult	7 ft O ₂ / 3 in. CO ₂	Female	25
1843	Softech® Bi-Flo ETCO ₂ Sampling Cannula, Adult	7 ft O ₂ / 3 in. CO ₂	Female	25
RHET44U	Softech® Bi-Flo ETCO ₂ Sampling Cannula, Adult	7 ft O ₂ / 7 ft CO ₂	Male	25
1844	Softech® Bi-Flo ETCO ₂ Sampling Cannula, Adult	7 ft O ₂ / 7 ft CO ₂	Male	25
RHET45U	Softech® Bi-Flo ETCO ₂ Sampling Cannula, Adult	7 ft O ₂ / 7 ft CO ₂	Female	25
1845	Softech® Bi-Flo ETCO ₂ Sampling Cannula, Adult	7 ft O ₂ / 7 ft CO ₂	Female	25

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.