


ETCO₂

SOFTECH® PLUS ETCO₂ SAMPLING CANNULA

Description:

The Softech® Plus ETCO₂ Sampling Cannula allows delivery of supplemental oxygen while simultaneously sampling expired CO₂, and features softer nasal prongs and lariat for enhanced patient comfort. The split nare design allows separation of the CO₂ readings and the delivery of oxygen, and helps produce a sharp wave form for clinician diagnostics.



United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699

Lat

Am

Catalog #	Size	Tubing Length	Luer Connector	Qty/Case
2842	Softech® Plus ETCO ₂ Sampling Cannula, Adult	7 ft O ₂ / 2 in. CO ₂	Male	50
2843	Softech® Plus ETCO ₂ Sampling Cannula, Adult	7 ft O ₂ / 2 in. CO ₂	Female	50
2844	Softech® Plus ETCO ₂ Sampling Cannula, Adult	7 ft O ₂ / 7 ft CO ₂	Female	50
2845	Softech® Plus ETCO ₂ Sampling Cannula, Adult	7 ft O ₂ / 7 ft CO ₂	Male	50
2846	Softech® Plus ETCO ₂ Sampling Cannula, Adult	14 ft O ₂ / 2 in. CO ₂	Female	50
2847	Softech® Plus ETCO ₂ Sampling Cannula, Adult	14 ft O ₂ / 14 ft CO ₂	Female	50

Catalog #	Size	Tubing Length	Luer Connector	Qty/Case
2848	Softech® Plus ETCO ₂ Sampling Cannula, Adult	14 ft O ₂ / 14 ft CO ₂	Male	50
2849	Softech® Plus ETCO ₂ Sampling Cannula, Pediatric	7 ft O ₂ / 2 in. CO ₂	Female	50
2850	Softech® Plus ETCO ₂ Sampling Cannula, Pediatric	7 ft O ₂ / 7 ft CO ₂	Female	50
2851	Softech® Plus ETCO ₂ Sampling Cannula, Pediatric	7 ft O ₂ / 7 ft CO ₂	Male	50
2852	Softech® Plus ETCO ₂ Sampling Cannula, Pediatric	14 ft O ₂ / 2 in. CO ₂	Female	50
2853	Softech® Plus ETCO ₂ Sampling Cannula, Pediatric	14 ft O ₂ / 14 ft CO ₂	Female	50
2854	Softech® Plus ETCO ₂ Sampling Cannula, Pediatric	14 ft O ₂ / 14 ft CO ₂	Male	50
2855	Softech® Plus ETCO ₂ Sampling Cannula, Infant	7 ft O ₂ / 2 in. CO ₂	Female	50
2856	Softech® Plus ETCO ₂ Sampling Cannula, Infant	7 ft O ₂ / 7 ft CO ₂	Female	50
2857	Softech® Plus ETCO ₂ Sampling Cannula, Infant	7 ft O ₂ / 7 ft CO ₂	Male	50

© 2026 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.