


GUEDEL ORAL AIRWAYS

TRADITIONAL GUEDEL AIRWAY

Description:

The Rusch® Traditional Guedel Airway features:

- Rigid plastic
- Economical choice
- Full range of sizes and convenient packaging options
- Size identified on the flange



United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699

Lat

Am

Catalog #	Size	Length	ISO Standard Size Design	Qty/Box
122250	0	50 mm	5	10 (Bulk)
122260	1	60 mm	6	10 (Bulk)
122270	2	70 mm	7	10 (Bulk)
122280	3	80 mm	8	10 (Bulk)
122290	4	90 mm	9	10 (Bulk)
1222100	5	100 mm	10	10 (Bulk)

Catalog #	Size	Length	ISO Standard Size Design	Qty/Box
122850	0	50 mm	5	50 (Bulk)
122860	1	60 mm	6	50 (Bulk)
122870	2	70 mm	7	50 (Bulk)
122880	3	80 mm	8	50 (Bulk)
122890	4	90 mm	9	50 (Bulk)
1228100	5	100 mm	10	50 (Bulk)
122150	0	50 mm	5	10 (Individually Packaged)
122160	1	60 mm	6	10 (Individually Packaged)
122170	2	70 mm	7	10 (Individually Packaged)
122180	3	80 mm	8	10 (Individually Packaged)
122190	4	90 mm	9	10 (Individually Packaged)
1221100	5	100 mm	10	10 (Individually Packaged)
122950	0	50 mm	5	50 (Individually Packaged)
122960	1	60 mm	6	50 (Individually Packaged)
122970	2	70 mm	7	50 (Individually Packaged)

Catalog #	Size	Length	ISO Standard Size Design	Qty/Box
122980	3	80 mm	8	50 (Individually Packaged)
122990	4	90 mm	9	50 (Individually Packaged)
1229100	5	100 mm	10	50 (Individually Packaged)

© 2026 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.