



Pilling®

## Maloney & Hurst-Style Esophageal Bougies


Silicone, Mercury-Free

### Maloney & Hurst-Style Esophageal Bougies

We offer two choices for esophageal dilation – tapered or blunt tip. Pilling Bougies are tungsten-filled and mercury-free, which means they can be steam sterilized and disposed of easily. The bougies are flexible but not made with natural rubber latex. The silicone casing is non-slip when dry

and slippery when wet. Pilling bougies give surgeons a good grip, but allow for easy insertion.


With a range of sizes to choose from, as well as a storage tray, contact your Teleflex Instrument Sales Specialist and let us help you with your esophageal bougie dilator needs.



### Maloney Tapered Tungsten-Filled Esophageal Bougies

This tungsten-filled bougie provides a single style of dilator for treating cardiospasm, esophagitis, stenosis, and other esophageal diseases. The small, flexible tip is easily introduced through the cricopharynx area, cardioesophageal junction or stenosis. Once the tip is through the involved area, the tapered end allows the single bougie to dilate this portion of the esophagus to its maximum limit. Overall length is approximately 29-1/2" (75 cm).

CATALOG #	SIZE	CATALOG#	SIZE	CATALOG#	SIZE
507920	20 Fr.	507934	34 Fr.	507948	48 Fr.
507922	22 Fr.	507936	36 Fr.	507950	50 Fr.
507924	24 Fr.	507938	38 Fr.	507952	52 Fr.
507926	26 Fr.	507940	40 Fr.	507954	54 Fr.
507928	28 Fr.	507942	42 Fr.	507956	56 Fr.
507930	30 Fr.	507944	44 Fr.	507958	58 Fr.
507932	32 Fr.	507946	46 Fr.	507960	60 Fr.



### Hurst Tungsten-Filled Esophageal Bougies

These tungsten-filled bougies have rounded distal tips. The closure on the proximal end is drilled through from side-to-side so a finger cord may be used. 29-1/2" (75 cm) approximate overall length.

CATALOG #	SIZE	CATALOG #	SIZE	CATALOG #	SIZE
507595	20 Fr.	507602	34 Fr.	507609	48 Fr.
507596	22 Fr.	507603	36 Fr.	507610	50 Fr.
507597	24 Fr.	507604	38 Fr.	507611	52 Fr.
507598	26 Fr.	507605	40 Fr.	507612	54 Fr.
507599	28 Fr.	507606	42 Fr.	507613	56 Fr.
507600	30 Fr.	507607	44 Fr.	507614	58 Fr.
507601	32 Fr.	507608	46 Fr.	507615	60 Fr.



### Bougie Tray

Made of Stainless steel. Separately holds 12 individual bougies (Maloney or Hurst) in two tiers - the upper tier with a hinged handle, perforated trays, and lid. Holding pins hold the bougie cap to keep bougie straight and in corrugated channels.

CATALOG #	DESCRIPTION
508000	31-5/8" x 11" (80.5 x 28 cm)



[teleflexsurgicalcatalog.com](http://teleflexsurgicalcatalog.com)

Teleflex, the Teleflex logo, and Pilling are trademarks or registered trademarks of Teleflex Incorporated or its affiliates, in the U.S. and/or other countries. Information in this document is not a substitute for the product Instructions for Use. The products in this catalogue may not be available in all countries. Please contact your local representative. All data current at time of printing (04/2017). Subject to technical changes without further notice. ©2017 Teleflex Incorporated. All rights reserved. MC-003020 Rev 0