



Pilling®
Low-Profile Dedo-Style Laser Laryngoscope
The Newest Addition to the Pilling® Family


Accommodating a Wider Range of Patient Anatomies

- Lower profile design
- Same internal diameter visualization
- Single light carrier and cable
- Matte finish inside, designed to be safer when using laser
- Same Pilling® quality



*Designed for use with Pilling Xenon 300W and XenaLight® XLT II
Lightsources from Teleflex (528300P & 528305XLTP, respectively)*

Compared to Previous-Generation Pilling® Dedo Scopes:



*522221 has an outer diameter that is 5% (1.5 mm) slimmer than 522225P and 20% (6.5 mm) slimmer than 522228

NEW Pilling® Low-Profile Laryngoscope


ITEM	DESCRIPTION	OVERALL LENGTH	MIN I.D.	NOMINAL PROXIMAL DIMENSIONS	MAX O.D.	REPLACEMENT LIGHT CARRIER
522221	Pilling® Low-Profile Dedo-Style Laser Laryngoscope	178 mm	10.7 mm	19.0 mm x 19.0 mm	25.5 mm	527868

Compared to Other Pilling® Dedo Designs

ITEM	DESCRIPTION	OVERALL LENGTH	MIN I.D.	NOMINAL PROXIMAL DIMENSIONS	MAX O.D.	REPLACEMENT LIGHT CARRIER
522225P	Pilling® Dedo Laryngoscope, Bifurcated Light Carrier	178 mm	10.7 mm	19.0 mm x 19.0 mm	27.0 mm	526525
522228	Pilling® Dedo Laser Laryngoscope, Dual Light Cables	177 mm	10.7 mm	20.0 mm x 19.0 mm	32.0 mm	527860
522231	Pilling® Ossoff-Karlan-Dedo Laser Laryngoscope, Dual Light Cables	178 mm	11.9 mm	22.0 mm x 20.0 mm	35.0 mm	527862

Pilling® ENT Double-Decker Instrument Storage Trays

ITEM	DESCRIPTION
P88097	ENT Double-Decker Tray with Mat and Inserts, 11" x 23" x 5"
P88044	ENT Double-Decker Tray with Inserts Only, 11" x 23" x 5"
087396	Ultra Soft Silicon Mat



teleflexsurgicalcatalog.com



Teleflex, the Teleflex logo, and Pilling are trademarks or registered trademarks of Teleflex Incorporated or its affiliates, in the U.S. and/or other countries. Information in this document is not a substitute for the product Instructions for Use. The products in this document may not be available in all countries. Please contact your local representative. All data current at time of printing (06/2019). Subject to technical changes without further notice.

© 2019 Teleflex Incorporated. All rights reserved. MC-005072 Rev 1

Teleflex®

Teleflex · 3015 Carrington Mill Boulevard, Morrisville, NC 27560
Toll Free: 866 246 6990 | Phone: +1 919 544 8000
teleflex.com