



**PILLING® SURGICAL LIGHTING SYSTEMS**  
Setting the Standard for Surgical Illumination Technology





# XENALIGHT™ XLT II LIGHTING SYSTEM

## XENALIGHT™ XLT II LIGHTSOURCE

- Five-port rotary color-coded turret
- Daylight-quality light
- Lamp lasts longer

## FIVE PORT ROTARY CONNECTION TURRET

- Color coded for ease of use
- Accepts the following endfittings: Wolf, Pilling®, Olympus, Storz, ACMI

## RIGIFLEX™ HEADBAND AND SURGILUX™ MODULE

- Lightweight and comfortable even in the longest procedures
- Available in fixed (100 mm) or variable (10-75 mm) spot

## FLOORSTAND

- Easy swivel; moves with the surgeon without resistance
- Rolling and locking casters
- Equipment hangers
- Handle that prevents equipment damage from accidental tip-over

## FIBEROPTIC CABLE

- Available in single or bifurcated




## DESIGNED FOR BRILLIANCE

For 200 Years, Pilling® Quality Instruments have been provided to surgeons across the world to improve the standard of patient care. Our portfolio of exceptional surgical lighting systems is part of that tradition.

In comparison to other illumination technologies, the XenaLight™ XLT II Lighting System offers the lowest decay curve over the life of the bulb. This equates to reduced light loss over the life of the bulb and can reduce maintenance time and costs over the lifetime of the lightsource.

**XENALIGHT™ XLT II LAMP GENERATES CONSISTENT WHITE LIGHT AND LONGER LAMP LIFE**



## XENALIGHT™ XLT II LIGHTSOURCE

- Specifically designed for Pilling® Laryngoscopes, Esophagoscopes and Bronchoscopes
- Designed to accommodate two Pilling Fiberoptic Cables

# LIGHTSOURCES SPECIFICALLY DESIGNED FOR YOUR PILLING® ENT SCOPES

## PILLING® LARYNGOSCOPIES



Offering one of the largest selections of laryngoscopes in the world today, Pilling Scopes include designs by leading surgeons. All laryngoscopes have fiberoptic light carriers.

## PILLING® TUCKER MEDIASTINOSCOPE



The improved Tucker Mediastinoscope has the highest light output of any mediastinoscope and is suited for use in conjunction with the operating microscope. The scope is supplied with an extra-large diameter Type II Fiberoptic light carrier and integral 6' long, autoclavable connecting cable.

## PILLING® PILLSBURY-JESBERG ESOPHAGOSCOPIES

The Pillsbury-Jesberg Esophagoscopes were designed with the light carrier fully incorporated into the textured handle. With this improved design, there is built-in protection for the light carrier during normal handling. The barrel and tip of the scope follow the familiar Jesberg design, facilitating gentle insertion while maintaining usual technique. Supplied with one light carrier. Available in multiple sizes. Blue fiberoptic cable sold separately.



## PILLING® IMPROVED HOLINGER VENTILATING FIBEROPTIC BRONCHOSCOPES

Our most popular bronchoscope, the Holinger, has been modified for pediatric foreign body removal and features the following:

- Small distal tip – the original Holinger tip has been cut down to correspond with tube size, while still retaining smoothness
- Built-in instrument channel for introduction of flexible forceps or suction
- Full lumen view – the oval fiberoptic light carrier has been raised slightly to provide a full lumen and give more room for instrumentation and a larger viewing area
- 15 mm sidearm attachment allows direct attachment of ventilation and anesthesia apparatus



524022	Telescope	2.7 mm, 0°
524036	Telescope Guide	2.7 mm
522552	Improved Holinger Bronchoscope	3.5 x 20
505670	Flexible 1 mm Alligator Forceps	

- Each Pilling Improved Holinger Bronchoscope is supplied with:
  - light carrier
  - one window plug
  - instrument guide rubber cap
  - blue fiberoptic cable sold separately

# FIBEROPTIC CABLES & ADAPTERS

## FIBEROPTIC CABLES

- Available in single and bifurcated cables
- 6' and 8' lengths
- Bifurcated cables allow simultaneous use of either two instruments or an instrument and headlight



## LIGHTSOURCE ENDFITTINGS



PILLING® LIGHTSOURCE ENDFITTING  
P356PRA, P358PRA



ACMI LIGHTSOURCE ENDFITTING  
A358PRA



OLYMPUS LIGHTSOURCE ENDFITTING  
O358PRA



STORZ LIGHTSOURCE ENDFITTING  
S358PRA



WOLF LIGHTSOURCE ENDFITTING  
W358PRA



BIFURCATED PILLING® BAYONET (BLUE SHEATHING)  
W356PXA, W358PXA

## INSTRUMENT ENDFITTINGS



ACMI FEMALE INSTRUMENT (AFA)



STORZ THREADED INSTRUMENT (STA)



PILLING® ROTATING INSTRUMENT (PRA)



**RIGIFLEX™ HEADBAND**  
 • Lightweight and comfortable  
 • Adjustable fit

**SURGILUX™ MODULE**  
 • Available in fixed (100 mm) or variable (10-75 mm) spot  
 • Adjustable focus

**FIBEROPTIC CABLE**  
 • Available in single or bifurcated

**POSITIONING CLIPS**  
 • Keeps cables securely fastened

**LIGHTSOURCE ENDFITTING**  
 • Wolf endfitting  
 • Compatible with a variety of adapters



**XENON 300W**

- Modern design that remains cool and runs quiet
- Patented five-port rotary color-coded turret
- Bulb hour meter and snap-in bulb cartridge offers easy maintenance

# ORDERING INFORMATION

LIGHTSOURCES	
PRODUCT #	DESCRIPTION
528300P	Xenon 300 Watt Lightsource, 5 Port Turret, 100-120V
528305XLTP	XenaLight™ XLT II Lightsource, 5 Port Turret, 100-120V


HEADLIGHT - SURGILUX™ MODULES	
PRODUCT #	DESCRIPTION
528341P	Variable 10-75 mm Spot with 8' XLG Cable
528342P	Fixed 100 mm Spot with 8' XLG Cable
528343P	Variable 10-75 mm Spot Headlight with 10' XLG Cable
528344P	Fixed 100 mm Spot with 10' XLG Cable
528339P	Variable 10-75 mm Spot with 8' XLG BIF Cable
528340P	Fixed 100 mm Spot with 8' XLG BIF Cable

FLOORSTANDS	
PRODUCT #	DESCRIPTION
528304P	Xenon 300 Watt, Base/Pole/Swivel
528309P	XenaLight™ XLT II Lighting System Base/Pole/Swivel

ACCESSORIES	
PRODUCT #	DESCRIPTION
528302P	Xenon 300 Watt Lightsource Bulb/Cartridge
528303P	Xenon 300 Watt Lightsource Bulb
528308P	XenaLight™ XLT Bulb with Timer Board for Older Model 528305, 528305XLT
528310P	XenaLight™ XLT II Bulb with Timer Chip for New Model 528305XLTP
W358XLGP	Headlight Cable 8', 3.5 mm, XLG Lensed (Wolf)
W3510XLGP	Headlight Cable 10', 3.5 mm, XLG Lensed (Wolf)
W358LWP	Headlight Cable 8', 3.5 mm, Bifurcated Lensed (Wolf)
W3510LWP	Headlight Cable 10', 3.5 mm, Bifurcated Lensed (Wolf)
528351P	Headband Front Insert, Breath-O-Prene (Pkg of 5)
528352P	Headband Top Insert, Breath-O-Prene (Pkg of 5)
528353P	Headband Back Insert, Breath-O-Prene (Pkg of 5)
528354P	Joysticks (Pkg of 5)
528355P	Gown Clips (Pkg of 10)



**ADAPTERS**

**LIGHTSOURCE END ADAPTER FOR FIBEROPTIC CABLE (FOR CUSTOMERS NOT USING A PILLING LIGHTSOURCE)**

50614AC	ACMI to Pilling® adapts to Pilling Single Cables with catalog numbers ending in "PRA" (formerly 521165)
50615AC	Wolf to Pilling® adapts to Pilling Single Cables with catalog numbers ending in "PRA" (formerly 521170)

**LIGHTSOURCE END ADAPTERS FOR TYPE II FIBEROPTIC CABLES - PILLING TO ACMI**

50616AC	<b>Adapts to the following:</b>
	527850 Light Carrier Cable for Improved Tucker Mediastinscope, long (522311)
	527860 Light Carrier Cable for Ossoff Adult Subglottiscopes (522244, 522245)
	527862 Light Carrier Cable for Ossoff-Karlan-Dedo Laser Micro-Laryngoscope, adol (522238)
	527863 Light Carrier Cable for Ossoff-Karlan-Dedo Laser Micro-Laryngoscope, child (522237)

**PILLING® ROTATING TIP LIGHT CARRIER CONNECTOR (IF MISPLACED)**

PRODUCT #	DESCRIPTION
50613AC	Adapts to Pilling® Single Cables with catalog numbers ending in "PRA" (formerly 521166)

FIBEROPTIC CABLES	INSTRUMENT ENDFITTINGS				
LIGHTSOURCE ENDFITTINGS	ACMI FEMALE INSTRUMENT	BIFURCATED (2) PILLING® ROTATING		PILLING® ROTATING BAYONET (PILLING BLUE)	STORZ THREADED INSTRUMENT
DESCRIPTION	8' LENGTH	6' LENGTH	8' LENGTH	6' LENGTH	8' LENGTH
Pilling® Lightsource Endfittings (P)				P356PRA	P358PRA
ACMI Lightsource Endfittings (A)	A358AFA				A358PRA
Olympus Lightsource Endfitting (O)					O358PRA
Storz Lightsource Endfitting (S)					S358PRA
Wolf Lightsource Endfittings (W)		W356PXA	W358PXA		W358PRA

Standard Bundle Size is 3.5 mm.  
 \* 5.0 mm Bundle Size

10' fiberoptic cables are available; contact your Instrument Sales Representative

- FIBEROPTIC CABLES**  
 Step 1: Select lightsource endfitting  
 Step 2: Select bundle size  
 Step 3: Select length  
 Step 4: Select instrument endfitting

Rigiflex, Surgilux and XenaLight are trademarks or registered trademarks of Aculux, Inc.

Teleflex and Pilling are trademarks or registered trademarks of Teleflex Incorporated or its affiliates.

Teleflex is a global provider of medical products designed to enable healthcare providers to protect against infections and improve patient and provider safety. The company specializes in products and services for vascular access, respiratory, general and regional anesthesia, cardiac care, urology and surgery. Teleflex also provides specialty products for device manufacturers.

© 2015 Teleflex Incorporated. All rights reserved. MC-000915



TELEFLEX 3015 Carrington Mill Boulevard, Morrisville, NC 27560  
Toll Free: 866.246.6990 Phone: +1.919.544.8000  
TELEFLEX.COM