

Catalog # :

502562

Albert-Sanders Venturi / Jet Ventilating Adapter for Laryngoscopes

Description:

Flow and pressure characteristics vary greatly according to the design of the bronchoscope or laryngoscope and any new combination should be studied carefully before clinical use.

Note: The Sanders Venturi System can be used only with tapered scopes.

Related:

502566 Sanders Venturi Adapter, 15 mm sidearm

NO IMAGE
AVAILABLE

United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699
