



**Bending Irons, Rod Holders, Ligature Carriers, Wire Passers****152255 : Deschamps Ligature Carrier**

Blunt point, left curved, 8-1/8" (20.5 cm) *Related: 152256 Blunt point, right curved, 8-1/8" (20.5 cm)*

**152256 : Deschamps Ligature Carrier**

Blunt point, right curved, 8-1/8" (20.5 cm) *Related: 152255 Blunt point, left curved, 8-1/8" (20.5 cm)*

**353640 : DeBakey® Modified Semb Ligature Carrier**

The distal 1.5 mm of the jaws are a DeBakey-type Design to securely hold a ligature or umbilical tape for passage around the aorta. Jaws curve on 2.5 cm radius. Curved, DeBakey Teeth, 10-3/8" (26.5 cm)

**353645 : Favalaro-Semb Ligature Carrier**

Jaws are constructed so that only the distal DeBakey- type Surfaces approximate. The smooth portion of the jaws remain open to prevent inadvertent clamping of the aortic wall. Half curved, DeBakey Teeth, 11-1/4" (28.5 cm)...

**353646 : Favalaro-Semb Ligature Carrier**

Jaws are constructed so that only the distal DeBakey- type Surfaces approximate. The smooth portion of the jaws remain open to prevent inadvertent clamping of the aortic wall. Full curved, DeBakey Teeth, 10-7/8" (27.5 cm)...

**KM46664 : UNIVERSAL ROD BENDER**

11 1/4" (28.5cm)



**KM46918 : Suture Passer with Suture Hole**

Used to pass sutures in shoulder surgery. Curved, 9-1/4" (23.5 cm). *Autoclavable to 290F (143C)* Related: KM46917 Suture Passer with Crochet Hook




**KM47919 : PLATE BENDING IRON**

for 3.5mm and 4.5mm plates




**KM48270 : KIRSCHNER BOW SMALL**

5 3/4" (14.5cm)




**KM48272 : KIRSCHNER BOW LARGE**

7 1/2" (19.0cm)



**KM48274 : KIRSCHNER BOW EXTRA LONG**

9 5/8" (24.5cm)




**KM48326 : WIRE AND PIN TRACTOR (BOEHLER) - SMALL W/2.4mm STEIN PIN**

Small, 90 mm inside length, 75 mm inside width. Pin size and diameter: 5" 2.4 mm. Includes pin KM71432. 6/pack.



**KM48328 : WIRE AND PIN TRACTOR (BOEHLER) - MEDIUM W/3.2mm STEIN PIN**

Medium, 160 mm inside length, 90 mm inside width. Pin size and diameter: 7" 3.2 mm. Includes pin KM71434. 6/pack.



**KM48330 : WIRE AND PIN TRACTOR (BOEHLER) - LARGE W/4.0mm STEIN PIN**

Large, 210 mm inside length, 150 mm inside width. Pin size and diameter: 9" 4.0 mm. Includes pin KM71525. 6/pack.




**KM48370 : GRATLOCH (TM) WIRE BENDER**

7 3/4" (19.5cm)




**KM48386 : Wire Passer Cannulated 45 mm**

Cannulation allows for easy feeding of wire around the bone. 45 mm diameter, 7" (18 cm). *Related: KM48388 70 mm diameter, 9-1/2" (24 cm)*




**KM48388 : Wire Passer Cannulated 70 mm**

Cannulation allows for easy feeding of wire around the bone. 70 mm diameter, 9-1/2" (24 cm). *Related: KM48386 45 mm diameter, 7" (18 cm)*



**KM49702 : ROD HOLDER 4.5MM THREADED**

8 1/4" (21.0cm)



**KM49706 : ROD HOLDER 5.0MM SMOOTH**

8 1/4" (21.0cm)



**KM49724 : K-WIRE & PLATE BENDER 3-PRONG MAX.1.1MM**

5" (12.5cm)



---

**KM49780 : IN SITU ROD BENDER 4.76MM RIGHT**

13" (33.0cm)



---

**KM49781 : IN SITU ROD BENDER 4.76MM LEFT**

13" (33.0cm)



---

**KM49782 : IN SITU ROD BENDER 6.35MM RIGHT**

13" (33.0cm)



---

**KM49783 : IN SITU ROD BENDER 6.35MM LEFT**

13" (33.0cm)



---

**KM49802 : K-WIRE BENDER 3 PRONG MAX.1.6MM**

5 3/4" (14.5cm)

