


K-Wire, Steinmann, Schanz Pins



KM71021 : KIRSCHNER WIRE

Smooth .028" (0.7mm) Double End Trocar, 6", BX/6



KM71022 : KIRSCHNER WIRE

Smooth .035" (0.9mm) Double End Trocar, 6", BX/6




KM71023 : KIRSCHNER WIRE

Smooth .045" (1.1mm) Double End Trocar, 6", BX/6




KM71024 : KIRSCHNER WIRE

Smooth .054" (1.4mm) Double End Trocar, 6", BX/6



KM71025 : KIRSCHNER WIRE

Smooth .062" (1.6mm) Double End Trocar, 6", BX/6



KM71031 : KIRSCHNER WIRE

Smooth .035" (0.9mm) Double End Diamond, 6", BX/6



KM71100 : KIRSCHNER WIRE

Smooth .028" (0.7mm) Double End Trocar, 4", BX/6



KM71101 : KIRSCHNER WIRE

Smooth .035" (0.9mm) Double End Trocar, 4", BX/6



KM71102 : KIRSCHNER WIRE

Smooth .045" (1.1mm) Double End Trocar, 4", BX/6



KM71103 : KIRSCHNER WIRE

Smooth .062" (1.4mm) Double End Trocar, 4", BX/6



KM71104 : KIRSCHNER WIRE

Smooth .054" (1.6mm) Double End Trocar, 4", BX/6



KM71106 : KIRSCHNER WIRE

Smooth .035" (0.9mm) Double End Trocar, 5", BX/6



KM71107 : KIRSCHNER WIRE

Smooth .045" (1.1mm) Double End Trocar, 5", BX/6





KM71108 : KIRSCHNER WIRE

Smooth .054" (1.4mm) Double End Trocar, 5", BX/6



KM71109 : KIRSCHNER WIRE

Smooth .062" (1.6mm) Double End Trocar, 5", BX/6



KM71110 : KIRSCHNER WIRE

Smooth .028" (0.7mm) Double End Trocar, 9", BX/6



KM71111 : KIRSCHNER WIRE

Smooth .035" (0.9mm) Double End Trocar, 9", BX/6



KM71112 : KIRSCHNER WIRE

Smooth .045" (1.1mm) Double End Trocar, 9", BX/6




KM71113 : KIRSCHNER WIRE

Smooth .062" (1.6mm) Double End Trocar, 9", BX/6




KM71114 : KIRSCHNER WIRE

Smooth .054" (1.4mm) Double End Trocar, 9", BX/6



KM71121 : KIRSCHNER WIRE

Smooth .035" (0.9mm) Double End Trocar, 12", BX/6



KM71131 : KIRSCHNER WIRE

Smooth .035" (0.9mm) Single Trocar, Round End, 4", BX/6



KM71132 : KIRSCHNER WIRE

Smooth .045" (1.1mm) Single Trocar, Round End, 4", BX/6



KM71133 : KIRSCHNER WIRE

Smooth .054" (1.4mm) Single Trocar, Round End, 4", BX/6



KM71134 : KIRSCHNER WIRE

Smooth .062" (1.6mm) Single Trocar, Round End, 4", BX/6

