


## Miscellaneous

**89-2718K : SP SAPH3 SS LIGHT PANEL, MINI****809319 : Weck® Disposable Electrodes**

Sterile, 24 per box. Blade electrode, 3/4" (19 mm) blade, 6" (15 cm) extended length.

**89-2717K : SP SAPH3 M LIGHT PANEL, MID****809337 : Weck® Disposable IMA/ENT Electrodes**


Sterile, 24 per box. Needle electrode, 3/4" (2 cm) needle, 1/32" diameter, extended insulation, 2-3/4" (7 cm).

**809317 : Weck® Disposable Electrodes**

Sterile, 24 per box. Needle electrode, 3/4" (19 mm) needle, 1/32" (0.7 mm) dia., 2-3/4" (7 cm).

**809335 : Weck® Disposable IMA/ENT Electrodes**

Sterile, 24 per box. Blade electrode, 3/4" (2 cm) blade, extended insulation, 2-3/4" (7 cm).



**89-2716K : SAPH III L LIGHT PANEL, LARGE**



**809339 : Weck® Disposable IMA/ENT Electrodes**

Sterile, 24 per box. Blade electrode, 3/4" (2 cm) blade, extended insulation, 6" (15 cm) extended length.



**351715 : CONN AORTIC COMPRESSOR**

15 1/2" (39.5cm)



**809318 : Weck® Disposable Electrodes**

Sterile, 24 per box. Needle electrode, 3/4" (19 mm) needle, 1/32" (0.7 mm) dia., 6" (15 cm) extended length.



**809338 : Weck® Disposable IMA/ENT Electrodes**

Sterile, 24 per box. Needle electrode, 3/4" (2 cm) needle, 1/32" diameter, extended insulation, 6" (15 cm) extended length.



**809316 : Weck® Disposable Electrodes**

Sterile, 24 per box. Needle electrode, 3/4" (19 mm) blunt needle, 1/32" (0.7 mm) dia., 2-3/4" (7 cm).



**809315 : Weck® Disposable Electrodes**

Sterile, 24 per box. Blade electrode, 3/4" (19 mm) blade, 2-3/4" (7 cm).



**522595 : WINDOW PLUG FOR BRONCHOSCOPES**

1 1/8" (2.5cm)



**601001A : UNIPOLAR BOVIE HIGH FREQUENCY CABLE**



**342002 : WOLVEK STERNAL APPROXIMATOR**

10 3/8" (26.5cm)




**507336 : LUER-LOCK CATHETER ADAPTER**

1 1/2" (4.0cm)



**524042 : INSTRUMENT GUIDE RUBBER CAPS**

5 per pack, 1/2" (1.5cm)



**522596 : REPLACEMENT LENS FOR WINDOW PLUG**

For 522595 window plug, 5/8" (1.5cm)



**429318 : UNIVERSAL EYE SHIELD**

20 per pack, 3 1/8" (8.0cm)



---

**342004 : WOLVEK STERNAL CRIMPER**



---

**507434 : ABRAHAM LARYNGEAL SYRINGE CANNULA W/LUER LOCK**

Malleable, 5 7/8" (15.0cm)



---

**502562 : Albert-Sanders Venturi / Jet Ventilating Adapter for Laryngoscopes**

Flow and pressure characteristics vary greatly according to the design of the bronchoscope or laryngoscope and any new combination should be studied carefully before clinical use. *Note: The Sanders Venturi System can be...*



---

**509214 : OPERATING PLATFORM**

Angled, Ebonized, 0.15" (4mm), 11 5/8" (29.5cm)



---

**509215 : OPERATING PLATFORM**

Angled, Ebonized, 0.25" (6mm), 11 5/8" (29.5cm)

