

Tunnelers, Tourniquets, Tendon Strippers, Retrievers**355497 : CRAWFORD-COOLEY GRAFT TUNNELER**


Graft is attached to obturator tip and drawn through the tube. Tube is then withdrawn leaving graft placed in untwisted position. Tube slightly curved, 12 mm I.D., 17-1/2" (44.5 cm)

**355500 : BERRY ANGIO ACCESS TUNNELER**


Designed by William B. Berry, M.D. of Pasadena, Texas, Angio Access Tunneler (355500) is used to create forearm loop tunnels for placement of arteriovenous grafts for hemodialysis. Removable interchangeable handle and...

**355495 : DEBAKEY® FEMORAL BYPASS TUNNELER**


For use in by-pass operation for segmental occlusive disease of the iliac and femoral arteries. This instrument is designed with a slightly malleable stainless steel bar with a handle on one end and an olive tip dilator on...

**065320 : KLEINERT-KUTZ TENDON RETRIEVER****065325 : CARROLL TENDON RETRIEVER**

Serrated

**353203 : STYLET FOR RUMEL TOURNIQUET**

10 1/4" (26cm)



353206 : STYLET FOR WECK STYLE RUMEL TOURNIQUET

9 3/8" (24cm)




353207 : STYLET CHILD SIZE RUMEL TOURNIQUET

7" (18cm)



KM46735 : KLEINERT-KUTZ TENDON RETRIEVER

Serrated, flexible shaft, 8" (20.3cm)



KM46739 : KLEINERT-KUTZ TENDON RETRIEVER

Serrated, 8" (20.3cm)



P18667 : CRAWFORD-COOLEY TUNNELER

12mm, 1.5" (4cm)