



**Forceps, Hemostats, Clamps****394236 : WECK® BIPOLAR FORCEPS CORDS**


Fits all standard bipolar forceps. Disposable, sterile, 12' (3.6 m), 10 per box.

**505100 : JACKSON LARYNGEAL CUP FORCEPS**


Straight, 4mm cup, 28 cm

**505102 : JACKSON LARYNGEAL CUP FORCEPS**

Straight, 6mm cup, 28 cm

**505104 : JACKSON LARYNGEAL CUP FORCEPS**


Angled, 4mm cup, 28 cm

**505106 : JACKSON LARYNGEAL CUP FORCEPS**

Angled, 6mm cup, 28 cm


**505230 : PATTERSON ELONGATED CUP BITING FORCEPS**

Straight, 11 1/8" (28 cm)




**506460 : JAKO-KLEINSASSER MICRO LARYNGEAL CUP FORCEPS**

Angled left, 2.4 mm tip, 8-5/8" (22 cm)



**841011 : KAISER-PILLING RING FORCEPS (2L)**

Long jaw, curved, 12-1/4" (31 cm)



**351802 : DEBAKEY TISSUE FORCEPS**

Straight, 2mm tip, 7 3/4" (20 cm)




**506472 : Jako-Kleinsasser Micro-Laryngeal Grasping Forceps**

These Pilling Micro-laryngeal Grasping Forceps were designed specifically for everyday-use micro-surgery. They feature slender 2 mm diameter stems with miniature distal jaws in the shape of grasping forceps and are made of...



**506455 : JAKO-KLEINSASSER MICRO LARYNGEAL CUP FORCEPS**

Straight, 2.4 mm tip, 8-5/8" (22 cm)



**506461 : JAKO-KLEINSASSER MICRO LARYNGEAL CUP FORCEPS**

Angled right, 2.4 mm tip, 8-5/8" (22 cm)



**505400 : JACKSON ALLIGATOR GRASPING FORCEPS**

Straight, Serrated, 23 x 3.8mm tip, 11 1/8" (28 cm)



**506462 : JAKO-KLEINSASSER MICRO LARYNGEAL CUP FORCEPS**

Angled up, 2.4 mm tip, 8-5/8" (22 cm)



**394246 : WECK® BICEPS FORCEPS CORDS**

Fits all Biceps Forceps. Disposable, sterile, 12' (3.6 m), 10 per box.




**443507 : JACOBSON HEMOSTATIC FORCEPS**

Curved, Serrated, 0.8mm tip, 5 1/8" (13 cm)




**211117A : LANDRENEAU MASHER GRASPING FORCEPS**

Lightly curved, 7.8 mm tip, 8-1/2" (22 cm)



**505116 : JACKSON LARYNGEAL CUP FORCEPS EXTRA LIGHT**

Angled, 2mm cup, 28 cm



**452550 : MAGILL CATHETER FORCEPS**

Serrated, 9 1/4" (24 cm)




**121655 : EDNA TOWEL FORCEPS**

Single row of non-penetrating teeth, 5-1/2" (14 cm)




**505112 : JACKSON LARYNGEAL CUP FORCEPS**

Straight, 2mm cup, 28 cm



**505420 : JACKSON ALLIGATOR GRASPING FORCEPS**

Straight, Serrated, 17 x 1.5mm tip, 11 3/4" (30 cm)



**841031 : KAISER-PILLING RING FORCEPS (2S)**

Curved 3-3/8" jaw, 1 x 11 ratchets, 12-1/8" (31 cm)



**182310 : HALSTEAD MOSQUITO FORCEPS**

Curved, Pilling-Pattern, 5 1/8" (13 cm)



**342981 : DEBAKEY THORACIC TISSUE FORCEPS**

Straight, 3mm tip, 7 3/4" (20 cm)

