



Forceps, Hemostats, Clamps

**394236 : WECK® BIPOLAR FORCEPS CORDS**


Fits all standard bipolar forceps. Disposable, sterile, 12' (3.6 m), 10 per box.

**505100 : JACKSON LARYNGEAL CUP FORCEPS**


Straight, 4mm cup, 28 cm

**505102 : JACKSON LARYNGEAL CUP FORCEPS**

Straight, 6mm cup, 28 cm

**505104 : JACKSON LARYNGEAL CUP FORCEPS**


Angled, 4mm cup, 28 cm

**505106 : JACKSON LARYNGEAL CUP FORCEPS**

Angled, 6mm cup, 28 cm


**505230 : PATTERSON ELONGATED CUP BITING FORCEPS**

Straight, 11 1/8" (28 cm)




506460 : JAKO-KLEINSASSER MICRO LARYNGEAL CUP FORCEPS

Angled left, 2.4 mm tip, 8-5/8" (22 cm)



841011 : KAISER-PILLING RING FORCEPS (2L)

Long jaw, curved, 12-1/4" (31 cm)



351802 : DEBAKEY TISSUE FORCEPS

Straight, 2mm tip, 7 3/4" (20 cm)




506472 : Jako-Kleinsasser Micro-Laryngeal Grasping Forceps

These Pilling Micro-laryngeal Grasping Forceps were designed specifically for everyday-use micro-surgery. They feature slender 2 mm diameter stems with miniature distal jaws in the shape of grasping forceps and are made of...




506455 : JAKO-KLEINSASSER MICRO LARYNGEAL CUP FORCEPS

Straight, 2.4 mm tip, 8-5/8" (22 cm)



506461 : JAKO-KLEINSASSER MICRO LARYNGEAL CUP FORCEPS

Angled right, 2.4 mm tip, 8-5/8" (22 cm)



505400 : JACKSON ALLIGATOR GRASPING FORCEPS

Straight, Serrated, 23 x 3.8mm tip, 11 1/8" (28 cm)



506462 : JAKO-KLEINSASSER MICRO LARYNGEAL CUP FORCEPS

Angled up, 2.4 mm tip, 8-5/8" (22 cm)



394246 : WECK® BICEPS FORCEPS CORDS

Fits all Biceps Forceps. Disposable, sterile, 12' (3.6 m), 10 per box.




443507 : JACOBSON HEMOSTATIC FORCEPS

Curved, Serrated, 0.8mm tip, 5 1/8" (13 cm)




211117A : LANDRENEAU MASHER GRASPING FORCEPS

Lightly curved, 7.8 mm tip, 8-1/2" (22 cm)



505116 : JACKSON LARYNGEAL CUP FORCEPS EXTRA LIGHT

Angled, 2mm cup, 28 cm



452550 : MAGILL CATHETER FORCEPS

Serrated, 9 1/4" (24 cm)




121655 : EDNA TOWEL FORCEPS

Single row of non-penetrating teeth, 5-1/2" (14 cm)




505112 : JACKSON LARYNGEAL CUP FORCEPS

Straight, 2mm cup, 28 cm



505420 : JACKSON ALLIGATOR GRASPING FORCEPS

Straight, Serrated, 17 x 1.5mm tip, 11 3/4" (30 cm)



841031 : KAISER-PILLING RING FORCEPS (2S)

Curved 3-3/8" jaw, 1 x 11 ratchets, 12-1/8" (31 cm)



182310 : HALSTEAD MOSQUITO FORCEPS

Curved, Pilling-Pattern, 5 1/8" (13 cm)



342981 : DEBAKEY THORACIC TISSUE FORCEPS

Straight, 3mm tip, 7 3/4" (20 cm)

