

Tracheal Laryngeal Tubes

**518035 : Jackson Improved Tracheostomy Tube #6**

Stainless steel. Interchangeable. The improved model with a rotating lock provides quick, easy removal of the inner tube. A set consists of an outer tube, inner tube and an obturator. 518035 (size 6, 10.0 mm O.D., 7.2 mm...

**518030 : Jackson Improved Tracheostomy Tube #5**

Stainless steel. Interchangeable. The improved model with a rotating lock provides quick, easy removal of the inner tube. A set consists of an outer tube, inner tube and an obturator. 518030 (size 5, 9.0 mm O.D., 6.2 mm...

**1037836 : DECANNULATION STOPPER**

White stopper with 40" cord, Size 6, Full Closure

**518025 : Jackson Improved Tracheostomy Tube #4**

Stainless steel. Interchangeable. The improved model with a rotating lock provides quick, easy removal of the inner tube. A set consists of an outer tube, inner tube and an obturator. 518025 (size 4, 8.0 mm O.D., 5.5 mm...

**518135 : JACKSON ORIGINAL TRACHEOSTOMY TUBE #6**

7.2mm I.D., 10mm O.D., 69mm

**518230 : JACKSON IMPROVED TRACHEOSTOMY TUBE #6, SHORT LENGTH**

7.2mm I.D., 10mm O.D., 56mm, stainless steel



518436 : JACKSON IMPROVED TRACHEOSTOMY TUBE #6, W/15MM ADAPTER, REGULAR LENGTH

7.2mm I.D., 10mm O.D., 69mm, stainless steel



518040 : Jackson Improved Tracheostomy Tube #7

Stainless steel. Interchangeable. The improved model with a rotating lock provides quick, easy removal of the inner tube. A set consists of an outer tube, inner tube and an obturator. 518040 (size 7, 11.0 mm O.D., 8.1 mm...



1037835 : DECANNULATION STOPPER

White stopper with 40" cord, Size 5, Full Closure



518902 : JACKSON TRACHEOSTOMY TUBE #6 W/ROTATING LOCK, EXTRA LONG

7.2mm I.D., 10mm O.D., 85mm, stainless steel



518225 : JACKSON IMPROVED TRACHEOSTOMY TUBE #5, SHORT LENGTH

6.2mm I.D., 9mm O.D., 56mm, stainless steel



1037834 : DECANNULATION STOPPER

White stopper with 40" cord, Size 4, Full Closure



518536 : JACKSON ORIGINAL TRACHEOSTOMY TUBE #6, W/15MM ADAPTER, REGULAR LENGTH

7.2mm I.D., 10mm O.D., 69mm, stainless steel



518431 : JACKSON IMPROVED TRACHEOSTOMY TUBE #5, W/15MM ADAPTER, REGULAR LENGTH

6.2mm I.D., 9mm O.D., 68mm, stainless steel



518045 : Jackson Improved Tracheostomy Tube #8

Stainless steel. Interchangeable. The improved model with a rotating lock provides quick, easy removal of the inner tube. A set consists of an outer tube, inner tube and an obturator. 518045 (size 8, 12.0 mm O.D., 9.2 mm...



518125 : JACKSON ORIGINAL TRACHEOSTOMY TUBE #4

5.5mm I.D., 8mm O.D., 62mm, stainless steel



518130 : JACKSON ORIGINAL TRACHEOSTOMY TUBE #5

6.2mm I.D., 9mm O.D., 68mm, stainless steel



518526 : JACKSON ORIGINAL TRACHEOSTOMY TUBE #4, W/15MM ADAPTER, REGULAR LENGTH

5.5mm I.D., 8mm O.D., 62mm, stainless steel



518220 : JACKSON IMPROVED TRACHEOSTOMY TUBE #4, SHORT LENGTH

5.5mm I.D., 8mm O.D., 52mm, stainless steel



518085 : SHIKANI-FRENCH SPEAKING VALVE FOR USE WITH JACKSON IMPROVED STAINLESS STEEL TRACHEOSTOMY TUBES

package of 2



**518250 : JACKSON IMPROVED TRACHEOSTOMY TUBE #10,
SHORT LENGTH**

10.6mm I.D., 14mm O.D., 56mm, stainless steel



**518900 : JACKSON TRACHEOSTOMY TUBE #5 W/ROTATING LOCK,
EXTRA LONG**

6.2mm I.D., 9mm O.D., 84mm, stainless steel



**518630 : JACKSON IMPROVED TRACHEOSTOMY TUBE #6,
W/15MM ADAPTER, SHORT LENGTH**

7.2mm I.D., 10mm O.D., 56mm, stainless steel



518039 : Jackson Improved Tracheostomy Tube - Fenestrated #6

Patients using standard Pilling Jackson Trachea Tubes can now realize the benefits provided by the Shikani-French Speaking Valve. By replacing the inner trachea tube with one of the tubes listed below, patients can attach...



**518426 : JACKSON IMPROVED TRACHEOSTOMY TUBE #4,
W/15MM ADAPTER, REGULAR LENGTH**

5.5mm I.D., 8mm O.D., 62mm, stainless steel