

Self-Retaining Retractors**341136 : PILLING® IMA RETRACTOR**

All components remain in the sterile field, giving you exposure and set up control.
Replaces 341134 Pilling Bugge IMA Retractor

**341139 : COUËTIL IMA STERNAL RETRACTOR****341195 : Cooley Sternum Retractor, X-Large**

Adult Spread 6-5/16" (16 cm) Long, concave blades 1-1/4" deep x 4" wide (3 x 10 cm) Arm length 6-1/2" (16.5 cm)

**351675 : Pilling® Multipurpose Retractor Complete**

A versatile self-retaining Henly-type Retractor designed for use in the following surgical specialties: • Peripheral vascular surgery in the carotid, subclavian, femoral, popliteal and tibial regions • Spinal surgery in...

**KM52771 : WILTSE-GELPI RETRACTOR**

Deep blades, sharp points, 10-7/8" (27.5 cm)

**341310 : Heitmiller Mini Thoracotomy Retractor**

Designed by Dr. Richard Heitmiller. 341310 This retractor provides excellent exposure for video-assisted thoracoscopy, open lung biopsies, transaxillary incisions and mid-CAB procedures in which mini-thoracotomy incisions...



341163 : Favaloro-Morse Sternal Retractor, Short

Spread 8" (20cm), Side Blade 1 1/8" (2.8cm) wide x 1 1/2" (3.8cm) deep, Arm length 6" (15.2cm)



341398 : KUYPER-MURPHY STERNAL RETRACTOR, MEDIUM

Concave blade, Spread 6 5/16" (16cm), Side Blade 3 1/2" (9cm) wide x 1 1/8" (3cm) deep, Straight arms



341180 : Gerbode Sternal Retractor

Spread 5 7/8" (15cm), Side Blade 3 15/16" (10cm) wide x 1 3/16" (3cm) deep



341312 : Premie Retractor

Spread 1 9/16" (4.0cm), Side Blade 1/4" (6mm) wide x 1/4" (6mm) deep, Arm length 1" (2.5cm)




341328 : FINOCHIETTO RETRACTOR, CHILD

Spread 5 1/8" (13cm), Side Blade 1 3/4" (4.4cm) wide x 1 1/2" (3.8cm) deep, Arm length 5 1/8" (13cm), Curved arms



351679 : IMPROVED HENLY RETRACTOR WITH STRAIGHT ARMS

7" (18cm), 3 pair of Side Blades 1" (2.5cm), 2" (5cm) and 3" (7.6cm) deep x 7/8" (2.2cm) wide. 3 Center Blades 7/8" (2.2cm), 1 1/8" (2.8cm), 1 1/4" (3.2cm) deep x 5/8" (1.6cm) wide, frame



382710 : ANDERSON ADSON RETRACTOR

Spread 3 1/2" (9cm), 4 x 4 Blunt prongs 17 X 7mm, 6 3/4" (17.1cm), Contoured arms



341150 : Chaux Multipurpose Sternal IMA Retractor

The Chaux Sternal Retractor is a unique instrument that can be used not only for routine midline sternotomy, but also for exposure and dissection of the left or right internal thoracic artery. In addition, accessory...



341317 : DeBakey® Child Retractor

Lightweight aluminum arms and blades are dull finished to eliminate glare and curved to follow body contour. Rack is stainless steel. Spread 3-3/4" (9.5 cm) Blades 1-1/8" deep x 1-1/8" wide (3 x 3 cm) Arm length 3-1/8" (8...



351680 : IMPROVED HENLY RETRACTOR WITH ARMS SPREAD OUTWARD

7" (18cm), 3 pair of Side Blades 1" (2.5cm), 2" (5cm) and 3" (7.6cm) deep x 7/8" (2.2cm) wide. 3 Center Blades 7/8" (2.2cm), 1 1/8" (2.8cm), 1 1/4" (3.2cm) deep x 5/8" (1.6cm) wide, frame



341194 : Cooley Sternum Retractor, Adult

Adult Spread 6-1/2" (16.5 cm) Short blades 1-1/4" deep x 2-1/4" wide (3.2 x 5.7 cm) Arm length 5-3/4" (14.6 cm)



341315 : Debakey® Infant Rib Spreader

A miniature model of the Finochietto. Light construction, stainless steel. Minimum spread 2-3/4" (7 cm), Blades 7/8" deep x 7/8" wide (22 mm deep x 21 mm wide)



485265 : MAHORNER THYROID RETRACTOR

7 Sharp prongs, 9" (23.0cm)



P17258 : ANKENY STERNAL RETRACTOR

6 sets of blades, 7" opening

AVAILABLE



341290 : Davis-Finochietto Retractor, Infant

Spread 2 3/4" (7cm), Side Blade 5/8" (1.6cm) wide x 1/2" (1.3cm) deep, Arm length 2 1/16" (5.2cm)

NO IMAGE AVAILABLE

341250 : Burford-Finochietto Child Retractor

Spread 5 7/8" (15cm), Side Blade 15/16" (2.4cm), 1 5/16" (3.3cm), 1 3/4" (4.4cm), 1 3/4" (4.4cm) wide x 1" (2.5cm), 1 1/4" (3.2cm), 1 1/4" (3.2cm), 1 15/16" (4.9cm) deep, Arm length 4 5/16" (11cm)



341302 : Finochietto Retractor, Infant

Spread 2 3/4" (7cm), Side Blade 9/16" (1.5cm) wide x 1/2" (1.3cm) deep, Arm length 2 1/8" (5.5cm), Aluminum



341181 : Gerbode Sternal Retractor

Spread 8 1/2" (21.5cm), Side Blade 3 15/16" (10cm) wide x 1 3/16" (3cm) deep



382722 : ADSON RETRACTOR

Spread 4 3/8" (11cm), 4 x 4 Sharp prongs 17 X 7mm, 7 1/4" (18.4cm), Contoured arms