

## Speculums

**212005 : GRAVES VAGINAL SPECULUM, MEDIUM**

Right open-sided vaginal specula, blade 4-3/8" long x 1-3/8" wide (11 x 3.5 cm)

**212185 : GARRIGUE WEIGHTED VAGINAL SPECULUM**

1 5/8" x 3 1/8" (4.0 cm x 8.0 cm), Weight 2 lbs.

**431120 : PECK-VIENNA NASAL SPECULUM**


5 3/8" (13.5cm)

**431111 : VIENNA NASAL SPECULUM, MEDIUM**

5 3/8" (13.5cm)

**212187 : HENROTIN WEIGHTED VAGINAL SPECULUM**

1 1/2" x 3 1/2" (3.8 cm x 9.0 cm), Weight 3 lbs.

**212011 : PEDERSON VAGINAL SPECULUM, MEDIUM**

4 3/8" (11.0 cm) x 1" (2.4 cm) blades



**212197 : AUVARD WEIGHTED VAGINAL SPECULUM**

Weight 3 lbs., 4" (10.2cm)



**212196 : AUVARD WEIGHTED VAGINAL SPECULUM**

Weight 2 1/2 lbs., 3 3/8" (8.6cm)



**431130 : COTTLE SEPTUM SPECULUM**

1 7/8" (5.0cm) blade, 5 1/2" (14.0cm)



**508614 : PILLING LASER SPECULUM, LARGE**

Blade 5" long x 1-1/2" wide (13 x 3.7 cm), suction tube 5-1/2" (14 cm) long (unbent), inner blades are matte finished to help reduce reflection of the laser beam



**030302 : KILLIAN NASAL SPECULUM, SIZE 2**

Blade 55mm, 5 1/8" (13.0cm)



**411140 : HARTMANN EAR SPECULUM**

Round, chrome plated, set of four



**025162 : FARRIER OVAL SPECULUM**

4 x 5mm



**508612 : PILLING LASER SPECULUM, MEDIUM**

Blade 4-3/8" long x 1-3/8" wide (11 x 3.5 cm), suction tube 4-5/8" (12 cm) long (unbent), inner blades are matte finished to help reduce reflection of the laser beam



**025165 : FARRIER OVAL SPECULUM**

5.5 x 6.5 mm



**030303 : KILLIAN NASAL SPECULUM, SIZE 3**

Blade 65mm, 5 1/8" (13.0cm)



**030304 : KILLIAN NASAL SPECULUM, SIZE 4**

Blade 80mm, 5 1/8" (13.0cm)



**202060 : PRATT RECTAL SPECULUM, LARGE**

Blades 3 3/8" (8.7 cm) long, 1" (2.5 cm) diameter



**030301 : KILLIAN NASAL SPECULUM, SIZE 1**

Blade 45mm, 5 1/8" (13.0cm)



**025164 : FARRIER OVAL SPECULUM**

5 x 6 mm



**411145 : HARTMANN EAR SPECULUM**

Oval, chrome plated, set of four



**508610 : PILLING LASER SPECULUM, SMALL**

Blade 3-1/4" long x 7/8" wide (8.5 x 2.2 cm), suction tube 4-5/8" (12 cm) long (unbent), inner blades are matte finished to help reduce reflection of the laser beam



**025166 : FARRIER OVAL SPECULUM**

6 x 7 mm



**025168 : FARRIER OVAL SPECULUM**

7 x 8mm



**202062 : PRATT RECTAL SPECULUM, SMALL**

Blades 3 1/2" (9.0cm) long, 3/4" (1.8cm) diameter