

Catheters

NO IMAGE
AVAILABLE**351470 : JEHLE CORONARY CATHETER, STERILE**

Cuff size diameter 9/32" (7.1mm), Catheter diameter .045" (1.14mm) I.D., .1" (2.54mm) O.D., overall length 20 1/4" (51.5cm)

NO IMAGE
AVAILABLE**351472 : JEHLE CORONARY CATHETER, STERILE**

Cuff size diameter 11/32" (8.7mm), Catheter diameter .045" (1.14mm) I.D., .1" (2.54mm) O.D., overall length 20 1/4" (51.5cm)

NO IMAGE
AVAILABLE**351474 : JEHLE CORONARY CATHETER, STERILE**

Cuff size diameter 13/32" (10.3mm), Catheter diameter .068" (1.73mm) I.D., .136" (3.45mm) O.D., overall length 20 1/4" (51.5cm)

NO IMAGE
AVAILABLE**351476 : JEHLE CORONARY CATHETER, STERILE**

Cuff size diameter 15/32" (11.9mm), Catheter diameter .068" (1.73mm) I.D., .136" (3.45mm) O.D., overall length 20 1/4" (51.5cm)

NO IMAGE
AVAILABLE**351490 : JEHLE CORONARY CATHETER, 4-HEAD, STERILE**

Cuff sizes 9/32", 11/32", 13/32", 15/32". Four tubes (2 tubes .1" O.D. x .045" I.D., 2 tubes .136" O.D. x .068" I.D.), overall length 20 1/2" (52cm)