

341195**341195 : Cooley Sternum Retractor, X-Large**

Adult Spread 6-5/16" (16 cm) Long, concave blades 1-1/4" deep x 4" wide (3 x 10 cm) Arm length 6-1/2" (16.5 cm)

**341506 : LEMOLE-PILLING® VALVE RETRACTOR - Set #1 with Cooley Retractor (341195)**

Provides self-retaining valve retraction. Fits directly onto arms of many types of sternal retractors. The T-Frame cross bar attaches to the sternal retractor by means of a "C" clamp. 8" (20.3 cm) T-Frame cross bar allows...