

Pilling® Blue Fiber Optic Cables**A358AFA : Pilling® Blue Fiber Optic Cable**

For use with Pilling Laryngoscopes, Bronchoscopes and Esophagoscopes. A358AFA: Pilling Blue Fiber Optic Cable, ACMI illuminator end, 3.5 mm bundle, 8 ft. length, ACMI Female instrument end

**A358PRA : Pilling® Blue Fiber Optic Cable**

For use with Pilling Laryngoscopes, Bronchoscopes and Esophagoscopes. A358PRA: Pilling Blue Fiber Optic Cable, ACMI illuminator end, 3.5 mm bundle, 8 ft. length, Pilling rotating instrument end (formally A358PR)

**O358PRA : Pilling® Blue Fiber Optic Cable**

For use with Pilling Laryngoscopes, Bronchoscopes and Esophagoscopes. O358PRA: Pilling Blue Fiber Optic Cable, Olympus illuminator end, 3.5 mm bundle, 8 ft. length, Pilling rotating instrument end

**P356PRA : Pilling® Blue Fiber Optic Cable**

For use with Pilling Laryngoscopes, Bronchoscopes and Esophagoscopes. P356PRA: Pilling Blue Fiber Optic Cable, Pilling illuminator end, 3.5 mm bundle, 6 ft. length, Pilling rotating instrument end

**P358PRA : Pilling® Blue Fiber Optic Cable**

For use with Pilling Laryngoscopes, Bronchoscopes and Esophagoscopes. P358PRA: Pilling Blue Fiber Optic Cable, Pilling illuminator end, 3.5 mm bundle, 8 ft. length, Pilling rotating instrument end

**S358PRA : Pilling® Blue Fiber Optic Cable**

For use with Pilling Laryngoscopes, Bronchoscopes and Esophagoscopes. S358PRA: Pilling Blue Fiber Optic Cable, Storz illuminator end, 3.5 mm bundle, 8 ft. length, Pilling rotating instrument end



W356PXA : Pilling® Bifurcated Endfitting Blue Fiberoptic Cable

For use with Pilling Laryngoscopes, Bronchoscopes and Esophagoscopes.W356PXA: Pilling Blue Fiber Optic Cable, Bifurcated cable, Wolf illuminator end, 3.5 mm bundle, 6 ft. length, Pilling rotating instrument end



W358PRA : Pilling® Blue Fiber Optic Cable

For use with Pilling Laryngoscopes, Bronchoscopes and Esophagoscopes.W358PRA: Pilling Blue Fiber Optic Cable, Wolf illuminator end, 3.5 mm bundle, 8 ft. length, Pilling rotating instrument end



W358PXA : Pilling® Bifurcated Endfitting Blue Fiberoptic Cable

For use with Pilling Laryngoscopes, Bronchoscopes and Esophagoscopes.W358PXA: Pilling Blue Fiber Optic Cable, Bifurcated cable, Wolf illuminator end, 3.5 mm bundle, 8 ft. length, Pilling rotating instrument end