

pilling right angle scope**341136 : PILLING® IMA RETRACTOR**

All components remain in the sterile field, giving you exposure and set up control.
Replaces 341134 Pilling Bugge IMA Retractor

**502265 : Dedo Table-Mounted Mayo Tray**


The Dedo Table-Mounted Mayo Tray System is used to support the standard Pilling Laryngoscope Holder (502248/502245). The tray can be positioned at any height, will swivel and can be extended. All parts are stainless steel...

**506476 : Jako-Kleinsasser Micro-Laryngeal Scissors**

These Pilling Micro-laryngeal Scissors were designed specifically for everyday-use micro-surgery. They feature a slender 2 mm diameter stem with miniature distal jaws in the shape of scissors and are made of dull finished...

**351675 : Pilling® Multipurpose Retractor Complete**

A versatile self-retaining Henly-type Retractor designed for use in the following surgical specialties: • Peripheral vascular surgery in the carotid, subclavian, femoral, popliteal and tibial regions • Spinal surgery in...

**506478 : Jako-Kleinsasser Micro-Laryngeal Scissors**

These Pilling Micro-laryngeal Scissors were designed specifically for everyday-use micro-surgery. They feature a slender 2 mm diameter stem with miniature distal jaws in the shape of scissors and are made of dull finished...

**354958 : Pilling® Thin Vascular and Coronary Technique® Needle Holder**

354958 Delicately proportioned, slender needle holders. 20.5 x 1.8 mm jaw, 9" (23 cm), TC



506475 : Jako-Kleinsasser Micro-Laryngeal Scissors

These Pilling Micro-laryngeal Scissors were designed specifically for everyday-use micro-surgery. They feature a slender 2 mm diameter stem with miniature distal jaws in the shape of scissors and are made of dull finished...



354962 : Pilling® Thin Vascular and Coronary Technique® Needle Holder

354962 Delicately proportioned, slender needle holders. 21.5 x 1.8 mm jaw, 8-1/4" (21 cm), TC




354961 : Pilling® Thin Vascular and Coronary Technique® Needle Holder

354961 Delicately proportioned, slender needle holders. 17 x 1.7 mm jaw, 7" (18 cm), TC



508614 : PILLING LASER SPECULUM, LARGE

Blade 5" long x 1-1/2" wide (13 x 3.7 cm), suction tube 5-1/2" (14 cm) long (unbent), inner blades are matte finished to help reduce reflection of the laser beam



354982 : PILLING EXTRA LIGHT NEEDLE HOLDER

14.5 x 1.6 mm jaw, 7" (18 cm), TC



506477 : Jako-Kleinsasser Micro-Laryngeal Scissors

These Pilling Micro-laryngeal Scissors were designed specifically for everyday-use micro-surgery. They feature a slender 2 mm diameter stem with miniature distal jaws in the shape of scissors and are made of dull finished...



508612 : PILLING LASER SPECULUM, MEDIUM

Blade 4-3/8" long x 1-3/8" wide (11 x 3.5 cm), suction tube 4-5/8" (12 cm) long (unbent), inner blades are matte finished to help reduce reflection of the laser beam



354959 : Pilling® Thin Vascular and Coronary Technique® Needle Holder

354959 Delicately proportioned, slender needle holders. 20.0 x 2.2 mm jaw, 10-3/8" (26.5 cm), TC



354960 : Pilling® Thin Vascular and Coronary Technique® Needle Holder

354960 Delicately proportioned, slender needle holders. 16.5 x 1.8 mm jaw, 6" (15 cm), TC



508610 : PILLING LASER SPECULUM, SMALL

Blade 3-1/4" long x 7/8" wide (8.5 x 2.2 cm), suction tube 4-5/8" (12 cm) long (unbent), inner blades are matte finished to help reduce reflection of the laser beam




518039 : Jackson Improved Tracheostomy Tube - Fenestrated #6

Patients using standard Pilling Jackson Trachea Tubes can now realize the benefits provided by the Shikani-French Speaking Valve. By replacing the inner trachea tube with one of the tubes listed below, patients can attach...




518034 : Jackson Improved Tracheostomy Tube - Fenestrated #5

Patients using standard Pilling Jackson Trachea Tubes can now realize the benefits provided by the Shikani-French Speaking Valve. By replacing the inner trachea tube with one of the tubes listed below, patients can attach...



352168 : Pilling Coronary Potts Scissors

Angled 120, 6-3/4" (17 cm), delicate, 16 mm blades




354981 : PILLING EXTRA LIGHT NEEDLE HOLDER

12.3 x 1.8mm jaw, 6" (15cm), TC



506479 : Jako-Kleinsasser Micro-Laryngeal Scissors

These Pilling Micro-laryngeal Scissors were designed specifically for everyday-use micro-surgery. They feature a slender 2 mm diameter stem with miniature distal jaws in the shape of scissors and are made of dull finished...



518049 : Jackson Improved Tracheostomy Tube - Fenestrated #8

Patients using standard Pilling Jackson Trachea Tubes can now realize the benefits provided by the Shikani-French Speaking Valve. By replacing the inner trachea tube with one of the tubes listed below, patients can attach...



351668 : RETRACTOR SIDE BLADE FOR 351675 PILLING MULTIPURPOSE RETRACTOR

4 Blunt prongs, Side Blade 7/8" (2.2cm) wide x 3 1/8" (8.0cm) deep



508004 : Video Laryngoscope Tray with Mat

Pilling Trays are designed for microsurgical and other delicate instruments. These lightweight, durable trays are engineered for maximum protection and flexibility. The transparent lid will help you easily identify...



352165 : Pilling Coronary Potts Scissors

Angled 45, 7" (18.0cm), Delicate, 16.0mm blades

