

retractor

NO IMAGE
AVAILABLE**211600 : Landreneau Thoroscopy Instrument Set**

211600: Complete set* includes all Landreneau catalog codes: 131114A, 211116A, 211117A, 211118A, 211603A, 211245A, 211290A, 211325A, 215558A, and 211500.

*Does not include 341310, the Heitmiller Mini Thoracotomy Retractor....

**341136 : PILLING® IMA RETRACTOR**

All components remain in the sterile field, giving you exposure and set up control.
Replaces 341134 Pilling Bugge IMA Retractor

**382805 : Dermahook "Fish Hook" Neuro Elastic Retractors**

1/2" (12.5 mm) hook, sterile, single use. 6/pouch, 10 pouches/box

**341139 : COUËTIL IMA STERNAL RETRACTOR****382800 : Durahook "Fish Hook" Neuro Elastic Retractors**

1/4" (6 mm) hook, sterile, single use. 6/pouch, 10 pouches/box

**341195 : Cooley Sternum Retractor, X-Large**

Adult Spread 6-5/16" (16 cm) Long, concave blades 1-1/4" deep x 4" wide (3 x 10 cm) Arm length 6-1/2" (16.5 cm)



341797 : Wylie Renal Vein Retractor

1" (2.5cm) width x 6 7/8" (17.5 cm) depth blade, 12 3/4" (32.5cm)



351675 : Pilling® Multipurpose Retractor Complete

A versatile self-retaining Henly-type Retractor designed for use in the following surgical specialties: • Peripheral vascular surgery in the carotid, subclavian, femoral, popliteal and tibial regions • Spinal surgery in...



343393 : Rail Clamps



KM52771 : WILTSE-GELPI RETRACTOR

Deep blades, sharp points, 10-7/8" (27.5 cm)



306010A : "S" Retractor

1/2" (1.5cm) width blade, 6 1/2" (16.5cm), 2/set



341310 : Heitmiller Mini Thoracotomy Retractor

Designed by Dr. Richard Heitmiller. 341310 This retractor provides excellent exposure for video-assisted thoracoscopy, open lung biopsies, transaxillary incisions and mid-CAB procedures in which mini-thoracotomy incisions...



341163 : Favaloro-Morse Sternal Retractor, Short

Spread 8" (20cm), Side Blade 1 1/8" (2.8cm) wide x 1 1/2" (3.8cm) deep, Arm length 6" (15.2cm)



381801A : O'Reilly Esophageal Retractor

3/8" (1.0cm) width blade, 21 1/4" (54.0cm)



341398 : KUYPER-MURPHY STERNAL RETRACTOR, MEDIUM

Concave blade, Spread 6 5/16" (16cm), Side Blade 3 1/2" (9cm) wide x 1 1/8" (3cm) deep, Straight arms




KM46838 : Hohmann Retractor

18mm width blade, 9" (23.0cm)



243200 : Lowsley Prostatic Tractor

Curved, 9 1/2" (24.0cm)



341180 : Gerbode Sternal Retractor

Spread 5 7/8" (15cm), Side Blade 3 15/16" (10cm) wide x 1 3/16" (3cm) deep




341312 : Premie Retractor

Spread 1 9/16" (4.0cm), Side Blade 1/4" (6mm) wide x 1/4" (6mm) deep, Arm length 1" (2.5cm)



342025 : Parsonnet Epicardal Retractor

3 X 3-pronged, 1 7/8" (4.5cm)





341328 : FINOCHIETTO RETRACTOR, CHILD

Spread 5 1/8" (13cm), Side Blade 1 3/4" (4.4cm) wide x 1 1/2" (3.8cm) deep, Arm length 5 1/8" (13cm), Curved arms



351779 : Henly Popliteal Retractor

4 X 3-pronged, Supplied w/3 center blades: 5/8" (1.5cm) width x 3/4" (1.8cm), 5/8" (1.5cm) x 1" (2.5cm) and 5/8" (1.5cm) x 1 1/4" (3.1cm) depth blade, 6 3/8" (16.0cm)



351679 : IMPROVED HENLY RETRACTOR WITH STRAIGHT ARMS

7" (18cm), 3 pair of Side Blades 1" (2.5cm), 2" (5cm) and 3" (7.6cm) deep x 7/8" (2.2cm) wide. 3 Center Blades 7/8" (2.2cm), 1 1/8" (2.8cm), 1 1/4" (3.2cm) deep x 5/8" (1.6cm) wide, frame



KM48220 : Saxena Bone Grafting Retractor

35mm graduated ratchet, 6 1/2" (16.5cm)



341795 : Wylie Renal Vein Retractor

1" (2.5cm) width x 3 3/4" (9.5cm) depth blade, 11 3/4" (30.0cm)