

Electrosurgery

**394236 : WECK® BIPOLAR FORCEPS CORDS**

Fits all standard bipolar forceps. Disposable, sterile, 12' (3.6 m), 10 per box.

**809319 : Weck® Disposable Electrodes**

Sterile, 24 per box. Blade electrode, 3/4" (19 mm) blade, 6" (15 cm) extended length.

**809337 : Weck® Disposable IMA/ENT Electrodes**

Sterile, 24 per box. Needle electrode, 3/4" (2 cm) needle, 1/32" diameter, extended insulation, 2-3/4" (7 cm).

**809317 : Weck® Disposable Electrodes**

Sterile, 24 per box. Needle electrode, 3/4" (19 mm) needle, 1/32" (0.7 mm) dia., 2-3/4" (7 cm).

**809335 : Weck® Disposable IMA/ENT Electrodes**

Sterile, 24 per box. Blade electrode, 3/4" (2 cm) blade, extended insulation, 2-3/4" (7 cm).

**809339 : Weck® Disposable IMA/ENT Electrodes**

Sterile, 24 per box. Blade electrode, 3/4" (2 cm) blade, extended insulation, 6" (15 cm) extended length.



809318 : Weck® Disposable Electrodes

Sterile, 24 per box. Needle electrode, 3/4" (19 mm) needle, 1/32" (0.7 mm) dia., 6" (15 cm) extended length.



809338 : Weck® Disposable IMA/ENT Electrodes

Sterile, 24 per box. Needle electrode, 3/4" (2 cm) needle, 1/32" diameter, extended insulation, 6" (15 cm) extended length.



809316 : Weck® Disposable Electrodes

Sterile, 24 per box. Needle electrode, 3/4" (19 mm) blunt needle, 1/32" (0.7 mm) dia., 2-3/4" (7 cm).




809315 : Weck® Disposable Electrodes

Sterile, 24 per box. Blade electrode, 3/4" (19 mm) blade, 2-3/4" (7 cm).




394246 : WECK® BICEPS FORCEPS CORDS

Fits all Biceps Forceps. Disposable, sterile, 12' (3.6 m), 10 per box.



392034 : WECK® BIPOLAR CUSHING BAYONET FORCEPS

Straight bipolar forceps, insulated, non-sterile. 0.7 mm tips, 7-1/2" (19 cm).



394258 : WECK® BICEPS BAYONET FORCEPS

Straight, non-stick bipolar forceps, non-sterile. 1.0 mm tips, with stops, 8" (20.5 cm).



392018 : WECK® BIPOLAR JEWELER'S FORCEPS

Straight, bipolar forceps, non-sterile. 0.4 mm tips, smooth handle, 4" (10 cm).



191025 : WECK® BIPOLAR CUSHING BAYONET FORCEPS

Straight bipolar forceps, insulated, non-sterile. 2.0 mm tips, 7-1/2" (19 cm).



392040 : WECK® BIPOLAR SEMKIN FORCEPS

Straight, bipolar forceps, non-sterile. 0.7 mm tips, with stops, 5-1/2" (14 cm).



392014 : WECK® BIPOLAR SCOVILLE-GREENWOOD BAYONET FORCEPS

Straight, bipolar forceps, non-sterile. 1.5 mm tips, insulated, 7-3/4" (19.5 cm).



392028 : WECK® BIPOLAR ADSON FORCEPS

Straight, bipolar forceps, non-sterile. 1.0 mm tips, insulated, 5" (12.5 cm).



392204 : WECK® BIPOLAR CUSHING BAYONET FORCEPS

Straight bipolar forceps, insulated, non-sterile. 1.5 mm tips, 7-3/4" (19.5 cm).



394811 : WECK® BICEPS CERULLO BAYONET FORCEPS

Straight, non-stick bipolar forceps, non-sterile. 0.5 mm tips, without stops, 10" (25.5 cm).





394238 : WECK® BIPOLAR FORCEPS CORDS

Fits all standard bipolar forceps. Silicone, non-sterile, 12' (3.6 m).



394237 : WECK® BIPOLAR FORCEPS CORDS

Fits all standard bipolar forceps. Rubber, non-sterile, 12' (3.6 m).



P42505 : ADSON BIPOLAR FORCEPS

Insulated, smooth jaws, non-stick



392026 : WECK® BIPOLAR ADSON FORCEPS

Straight, bipolar forceps, non-sterile. 1.0 mm tips, non-insulated, 5" (12.5 cm).



394807 : WECK® BICEPS BAYONET FORCEPS

Straight, non-stick bipolar forceps, non-sterile. 0.5 mm tips, without stops, 7-1/2" (19 cm).

