

**507952 : Maloney Tapered Tungsten-Filled Esophageal Bougie**

This tungsten-filled bougie provides a single style of dilator for treating cardiospasm, esophagitis, stenosis and other esophageal diseases. The small, flexible tip is easily introduced through the cricopharynx area,...

**508000 : Bougie Storage Tray - Holds 12 Bougies**


Stainless steel. Separately holds twelve (12) individual bougies (Maloney or Hurst). Two tiers with upper tier having hinged handle, perforated trays and lid. Holding pins hold bougie cap to keep bougie straight and in...

**518035 : Jackson Improved Tracheostomy Tube #6**

Stainless steel. Interchangeable. The improved model with a rotating lock provides quick, easy removal of the inner tube. A set consists of an outer tube, inner tube and an obturator. 518035 (size 6, 10.0 mm O.D., 7.2 mm...

**502265 : Dedo Table-Mounted Mayo Tray**

The Dedo Table-Mounted Mayo Tray System is used to support the standard Pilling Laryngoscope Holder (502248/502245). The tray can be positioned at any height, will swivel and can be extended. All parts are stainless steel...

**505100 : JACKSON LARYNGEAL CUP FORCEPS**

Straight, 4mm cup, 28 cm

**507956 : Maloney Tapered Tungsten-Filled Esophageal Bougie**

This tungsten-filled bougie provides a single style of dilator for treating cardiospasm, esophagitis, stenosis and other esophageal diseases. The small, flexible tip is easily introduced through the cricopharynx area,...



518030 : Jackson Improved Tracheostomy Tube #5

Stainless steel. Interchangeable. The improved model with a rotating lock provides quick, easy removal of the inner tube. A set consists of an outer tube, inner tube and an obturator. 518030 (size 5, 9.0 mm O.D., 6.2 mm...




1037836 : DECANNULATION STOPPER

White stopper with 40" cord, Size 6, Full Closure



431120 : PECK-VIENNA NASAL SPECULUM

5 3/8" (13.5cm)



505102 : JACKSON LARYNGEAL CUP FORCEPS

Straight, 6mm cup, 28 cm



506476 : Jako-Kleinsasser Micro-Laryngeal Scissors

These Pilling Micro-laryngeal Scissors were designed specifically for everyday-use micro-surgery. They feature a slender 2 mm diameter stem with miniature distal jaws in the shape of scissors and are made of dull finished...




507960 : Maloney Tapered Tungsten-Filled Esophageal Bougie

This tungsten-filled bougie provides a single stylet of dilator for treating cardiospasm, esophagitis, stenosis and other esophageal diseases. The small, flexible tip is easily introduced through the cricopharynx area,...



431111 : VIENNA NASAL SPECULUM, MEDIUM

5 3/8" (13.5cm)



505104 : JACKSON LARYNGEAL CUP FORCEPS

Angled, 4mm cup, 28 cm




507950 : Maloney Tapered Tungsten-Filled Esophageal Bougie

This tungsten-filled bougie provides a single style of dilator for treating cardiospasm, esophagitis, stenosis and other esophageal diseases. The small, flexible tip is easily introduced through the cricopharynx area,...



518025 : Jackson Improved Tracheostomy Tube #4

Stainless steel. Interchangeable. The improved model with a rotating lock provides quick, easy removal of the inner tube. A set consists of an outer tube, inner tube and an obturator. 518025 (size 4, 8.0 mm O.D., 5.5 mm...



505106 : JACKSON LARYNGEAL CUP FORCEPS

Angled, 6mm cup, 28 cm



506478 : Jako-Kleinsasser Micro-Laryngeal Scissors

These Pilling Micro-laryngeal Scissors were designed specifically for everyday-use micro-surgery. They feature a slender 2 mm diameter stem with miniature distal jaws in the shape of scissors and are made of dull finished...




507954 : Maloney Tapered Tungsten-Filled Esophageal Bougie

This tungsten-filled bougie provides a single style of dilator for treating cardiospasm, esophagitis, stenosis and other esophageal diseases. The small, flexible tip is easily introduced through the cricopharynx area,...




518135 : JACKSON ORIGINAL TRACHEOSTOMY TUBE #6

7.2mm I.D., 10mm O.D., 69mm



505230 : PATTERSON ELONGATED CUP BITING FORCEPS

Straight, 11 1/8" (28 cm)



506475 : Jako-Kleinsasser Micro-Laryngeal Scissors

These Pilling Micro-laryngeal Scissors were designed specifically for everyday-use micro-surgery. They feature a slender 2 mm diameter stem with miniature distal jaws in the shape of scissors and are made of dull finished...



507948 : Maloney Tapered Tungsten-Filled Esophageal Bougie

This tungsten-filled bougie provides a single style of dilator for treating cardiospasm, esophagitis, stenosis and other esophageal diseases. The small, flexible tip is easily introduced through the cricopharynx area,...



514413 : JACKSON TENACULUM TRACHEAL

Hook depth 1/4" (0.5 cm), 5-1/4" (13.5 cm)



518230 : JACKSON IMPROVED TRACHEOSTOMY TUBE #6, SHORT LENGTH

7.2mm I.D., 10mm O.D., 56mm, stainless steel

