

ACCESS

Balloon Access Ports

Description:

Teleflex has a broad offering in balloon ports predominantly used on primary entry utilizing open technique. Weck Vista® balloon access ports are available in five options all of which include a foam collar with locking clamp for easy one hand positioning. The foam collar is designed to maintain pneumoperitoneum. *Products listed contain no natural rubber latex.*



Catalog #	Details	Qty/Box
410944	WECK Vista® BALLOON ACCESS PORT 5/10 mm x 70 mm	6/BX
410944L	WECK Vista® BALLOON ACCESS PORT 5/10 mm x 100 mm	6/BX
410944S	WECK Vista® BALLOON ACCESS PORT 5/10 mm x 53 mm	6/BX
412944	WECK Vista® BALLOON ACCESS PORT 5/10/12 mm x 70 mm	6/BX
412944L	WECK Vista® BALLOON ACCESS PORT 5/10/12 mm x 100 mm	6/BX