


CHOLANGIOGRAPHY

Common Bile Duct Exploration

Description:

The Taut® System ONE™ CBDE Kit is an integrated combination of products designed specifically for transcystic laparoscopic common bile duct exploration. Tools needed to extract or flush CBD stones are provided in one convenient kit. Disposable components within each kit may be ordered. *Products listed contain no natural rubber latex.*



United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699

Lat

Am

Catalog #	Details	Qty/Box
50000	TAUT® SYSTEM ONE™ COMPREHENSIVE LAP CBDE KIT Includes 50065, 50640, 50024, PI-128, 50035 and K05-01251	1/BX
50024	TAUT® SYSTEM ONE™ TITANIUM4™ RETRIEVAL BASKET 4 flat wire titanium alloy basket with 360 degree articulation. Wire 2.4 French / 0.8 mm diameter. Open basket 15 mm dia. Overall length 100 cm	1/BX
50035	TAUT® SYSTEM ONE™ SOFTWIRE™ GUIDEWIRE Dual flexible ended PTFE coated guidewire. Wire 0.035 in. dia. Overall length 125 cm / 49.2 in.	1/BX
50065	TAUT® SYSTEM ONE™ REUSABLE ORGANIZER	1/BX
50640	TAUT® SYSTEM ONE™ MULTIPURPOSE BALLOON CATHETER Balloon catheter with radiopaque markers for clear visualization of distal and proximal balloon ends. Catheter 5 French / 1.67 mm dia., inflated 6 mm dia. Balloon length 4 cm. Overall length 75 cm	1/BX
K05-01251	TAUT® SYSTEM ONE™ HIGH PRESSURE INFLATION SYRINGE High pressure syringe for inflation of balloon catheter	1/BX