


**ACCESS****Cone Open Access Ports (Hasson Type)****Description:**

For those surgeons whose preference for primary entry is an open technique: Teleflex offers 12 mm diameter Weck Vista® ports. A locking clamp is included as part of the cone which allows single handed positioning of the cone. *Products listed contain no natural rubber latex.*



Catalog #	Details	Qty/Box
405933	WECK Vista® CONE OPEN ACCESS (HASSON TYPE) 5/10/12 mm x 100 mm	6/BX
405933L	WECK Vista® CONE OPEN ACCESS (HASSON TYPE) 5/10/12 mm x 125 mm	6/BX