

**LIGATION SYSTEMS**

**Horizon™ Metal Ligation System - Clips**

**Description:**

The preformed chevron shape of the Horizon™ Clip engulfs tissue with precise tip-to-tip closure. Heartshaped Titanium wire is designed to give each clip a firm grip on vessels, while the triangulated cross-section of the clip leg maximizes surface-to-surface contact between clip and jaw, practically eliminating clip fallout.



United States

To order this product, please call:

**866-246-6990**











Canada


To order this product, please call:

**800-387-9699**

Lat

Am

Catalog #	Clip Size	Clips/Cartridge	Cartridges/Box	Clips/Box
 005200	Microclip	6	30	180
 001200	Small	6	30	180
 001204	Small	24	25	600
 001201	Small-Wide	6	30	180
 001205	Small-Wide	24	25	600
 002200	Medium	6	30	180
 002204	Medium	24	25	600
 003200	Medium-Large	6	20	120
 003204	Medium-Large	24	15	360
 004200	Large	6	20	120

<b>Catalog #</b>	<b>Clip Size</b>	<b>Clips/Cartridge</b>	<b>Cartridges/Box</b>	<b>Clips/Box</b>
 004204	Large	24	15	360

© 2022 Teleflex Incorporated. All rights reserved.