

LIGATION SYSTEMS

Weck® Hem-o-lok® Polymer Ligation System - Appliers

Description:

The Weck® Hem-o-lok® Closure System's unique design and flexible hinge keep the clip firmly seated in applier jaws. The clip's distinct tactile feedback confirms jaw seating and secure vessel placement. *Weck® Hem-o-lok® appliers have been designed for exclusive use with Weck® Hem-o-lok® ligating clips.*



United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699

Lat

Am

Weck® Hem-o-lok® Polymer Open Appliers and Removers

Weck® Hem-o-lok® Appliers are available with curved or right-angled (70°) jaws and a wide range of applier lengths. Appliers open at box-lock for easy, thorough cleaning.

Catalog #	Clip Size	Applier Length	Jaw Type
■ 544113	Medium	8" (20 cm)	Curved
■ 544114	Medium	11" (28 cm)	Right-Angled (70°)
■ 544115	Medium	11" (28 cm)	Curved
■ 544123	Medium, Medium-Large, Large	11" (28 cm)	Remover
■ 544170	Medium-Large	8" (20 cm)	Curved
■ 544171	Medium-Large	11" (28 cm)	Curved

Catalog #	Clip Size	Applier Length	Jaw Type
 544172	Medium-Large	11" (28 cm)	Right-Angled (70°)
 544123	Medium-Large, Medium, Large	11" (28 cm)	Remover
 544179	Large	11" (28 cm)	Right-Angled (70°)
 544180	Large	8" (20 cm)	Curved
 544181	Large	11" (28 cm)	Curved
 544182	Large	10" (25.5 cm)	Right-Angled (70°)
 544123	Large, Medium, Medium-Large	11" (28 cm)	Remover
 544124	Large, Extra-Large	11" (28 cm)	Remover
 544191	Extra-Large	11" (28 cm)	Curved
 544192	Extra-Large	11" (28 cm)	Right-Angled (70°)
 544124	Extra-Large, Large	11" (28 cm)	Remover