


## LIGATION SYSTEMS

### Weck® Hem-o-lok® Polymer Ligation System - Clips

#### Description:

Weck® Hem-o-lok® Non-absorbable Polymer Locking Clips provide cool secure ligation; clip locks onto the vessel. The Weck® Hem-o-lok® system is available in a range of clip sizes Medium, Medium-Large, Large and Extra-Large. The Non-Absorbable polymer clip is inert, nonconductive, radiolucent, and does not interfere with CT or X-ray diagnostics.



United States

To order this product, please call:





**866-246-6990**

Canada

To order this product, please call:

**800-387-9699**

Lat  
Am

Catalog #	Clip Size	Clips/Cartridge	Cartridges/Box	Clips/Box
 544220	Medium	6	14	84
 544230	Medium-Large	6	14	84
 544240	Large	6	14	84
 544250	Extra-Large	6	14	84