

Home // Airway Management // Endotracheal Tubes // Specialty Intubation // Rusch® LaserTubus™ Magill Cuffed ETT  
// 102004-000070

CATALOG # : 102004-000070

## RUSCH® LASER TUBE LASER RESISTANT ETT

### Description:

Rusch® Laser Tube Laser Resistant Magill Cuffed ETT, Size 7.0, 2/BX

NO IMAGE  
AVAILABLE

United States

To order this product, please call:

**866-246-6990**

Canada

To order this product, please call:

**800-387-9699**

© 2026 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.