

Home // Airway Management // Endotracheal Tubes // Specialty Intubation // Rusch® Reinforced Uncuffed ETT // 103600-000040

CATALOG # : 103600-000040

RUSCH® REINFORCED UNCUFFED ETT

Description:

Rusch® Reinforced Uncuffed ETT, Size 4.0, 2/BX

NO IMAGE
AVAILABLE

United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699

© 2026 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.