


Home // Airway Management // Endotracheal Tubes // Standard Intubation // Rusch® Super Safety™ Magill Cuffed ETT // 112480035

CATALOG # : 112480035

## RUSCH® SUPER SAFETY CLEAR CUFFED ETT

### Description:

Rusch® Super Safety Clear Magill Cuffed ETT, Size 3.5, 10/BX



United States

To order this product, please call:

**866-246-6990**

Canada

To order this product, please call:

**800-387-9699**

Lat

Am

© 2026 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.