


Home // Airway Management // Endotracheal Tubes // Endotracheal Tube Accessories // Rusch® EndoTest™ Cuff Pressure Manometer // 112700-000000

CATALOG # : 112700-000000

RUSCH® ENDOTEST™ CUFF PRESSURE MANOMETER

Description:

Rusch® EndoTest™ Cuff Pressure Manometer, 1/BX



United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699

Lat

Am

© 2026 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.