


Home // Airway Management // Airways // Adjustable Flange Nasal Airways // Adjustable Flange Red Rubber Airway
// 125200-000240

CATALOG # : 125200-000240

ADJUSTABLE FLANGE RED RUBBER AIRWAY

Description:

Rusch® Adjustable Flange Red Rubber Airway 24 Fr x 170 mm, 5/BX (Bulk, non-sterile)



NO IMAGE
AVAILABLE

United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699

© 2026 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.