


Home // Airway Management // Airways // Adjustable Flange Nasal Airways // Adjustable Flange Red Rubber Airway
// 125200-000340

CATALOG # : 125200-000340

ADJUSTABLE FLANGE RED RUBBER AIRWAY

Description:

Rusch® Adjustable Flange Red Rubber Airway 34 Fr x 170 mm, 5/BX (Bulk, non-sterile)



United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699

© 2026 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.