

Home // Airway Management // LMA® Airways // LMA® Airway Accessories // LMA® Fastrach® Reusable Airway Accessories // 130099

CATALOG # : 130099

LMA FASTRACH® REUSABLE AIRWAY ACCESSORIES

Description:

LMA Fastrach® Reusable Stabilizer Rod, All ETT Sizes, 1/PK



United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699

© 2026 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.