

Home // Airway Management // LMA® Airways // Reusable LMA® Airways // LMA® Fastrach® Airway // 131065

CATALOG # : 131065

LMA FASTRACH® REUSABLE ETT

Description:

LMA Fastrach® Reusable ETT, Size 6.5 mm, 1/PK

NO IMAGE
AVAILABLE

United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699

Details:

- LMA Fastrach® Reusable ETT and Stabilizer Rod

© 2026 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.