

Home // Airway Management // Laryngoscopes // Reusable Laryngoscopes // Rüsch® Reusable Standard Handles  
// 142600

CATALOG # : 142600

## RUSCH® MATTE FINISH HANDLES

### Description:

Rusch® Matte Finish Standard Reusable Laryngoscope Handle, Small/Penlight, 1/BX



United States

To order this product, please call:

**866-246-6990**

Lat

Am

© 2026 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.