

Home // Airway Management // Humidification // Heat & Moisture Exchangers (HME) // Gibeck® Humid-Vent® HME
// 14412

CATALOG # : 14412

HUMID-VENT 2S

Description:

Humid-Vent 2S



To order this product, please call:
+1 (919) 433-4999

© 2026 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.