

Home // Airway Management // Laryngoscopes // Reusable Laryngoscopes // RUSCH® SnapLight® Reusable Fiber Optic Blades // 2515000

CATALOG # : 2515000

RUSCH® SNAPLIGHT™ BLADE

Description:

Rusch® SnapLight™ Fiber Optic Replacement Bundle, fits Item Code 004458800 (Miller OxiPort 0), 2/BX



United States

To order this product, please call:

866-246-6990

Lat

Am

© 2026 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.