


Home // Airway Management // Laryngoscopes // Reusable Laryngoscopes // RUSCH® SnapLight® Reusable Fiber Optic Blades // 2518300

CATALOG # : 2518300

## RUSCH® SNAPLIGHT™ BLADE

### Description:

Rusch® SnapLight™ Fiber Optic Replacement Bundle, fits Item Code 002218300 (Miller 3), 2/BX



United States

To order this product, please call:

**866-246-6990**

Lat  
Am

© 2026 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.