

CATALOG # : 2852

RESPIRATORY

Description:

Softech® Plus ETCO₂ Sampling Cannula, Pediatric, 50/CS (Individually Packaged)



United States

To order this product, please call:
866-246-6990

Lat
Am

- Details:
- Pediatric Cannula with 14 ft oxygen supply tubing and 2 in. CO₂ sampling line with a female luer connector

© 2026 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.